
anesthPlot

Release beta

Christophe Desbois

Jan 13, 2022

MAIN_sCRIP*T* :

1	Welcome to anesthPlot's documentation!	1
2	main script	3
3	modules	5
4	Indices and tables	21
	Python Module Index	23
	Index	25

WELCOME TO ANESTHPLOT'S DOCUMENTATION!

anesthPlot is a python package developped to extract, manipulate and plots anesthesia data recorded from the Monitor Software to be used mostly in a teaching environment.

Warning: This project is:

- a work in progres
- the processes are mainly focused on horses anesthesia
- in our environment the data recorded came from an as3 or as5 anesthesia machine monitoring ekg, invasive pressure, etCO2, halogenate, spirometry.

1.1 Features

- **load** recordings from a trend or a wave recordings

– **from command line:**

```
python anesthPlot/anesplot/__main__.py
```

- * build a **standard debriefing** (trends) **plot series** (script usage)
 - global histograms (cardiovascular and anesthesia summary)
 - cardiovascular trends time based plots
 - respiratory trends time based plots
 - anesthesia trends time based plots
- * build a **plot for wave** recording
 - one or two waves on the same plot (script usage)

- can also be used as a **python package**

– **usage :**

```
import anesplot.record_main as rec
trendname = 'a_full_path_to_csv_file'
trends = rec.MonitorTrend(trendname)
wavename = rec.trendname_to_wavename(trendname)
waves = rec.MonitorWave(trends)
```

(continues on next page)

(continued from previous page)

```
trends.show_graphs() # -> set of plots for debriefing purposes  
waves.plot_waves() # -> one or two traces  
waves.define_a_roi() # -> to register the plotting scales  
waves.animate_fig() #-> to build an animation using these parameters
```

- additional functions are available to extract instantaneous heart rate
 - * see `anesplot/treatrec/ekg_to_hr.py`

MAIN SCRIPT

2.1 anesplot.record_main module

main script/module to load and display an anesthesia record

can be runned as a script:: python record_main.py

or imported as a package:: import anesplot.record_main as rec %gui qt5 (required only to use the dialogs if using spyder) trends = rec.MonitorTrend() waves = rec.MonitorWave(rec.trendname_to_wavename(trends.filename))

anesplot.record_main.**choosefile_gui**(dirname=None)

Select a file via a dialog and return the (full) filename.

Parameters **dir_path** (str) – location to place the gui ('generally paths['data']) else home

Returns **fname[0]** – filename

Return type str

anesplot.record_main.**trendname_to_wavename**(name)

just compute the supposed name

anesplot.record_main.**select_type**(question=None, items=None, num=0)

select the recording type:

Returns **kind** – kind of recording in [monitorTrend, monitorWave, taphTrend, telvet]

Return type str

anesplot.record_main.**select_wave**(waves, num=1)

select the recording type:

Returns **kind** – kind of recording in [monitorTrend, monitorWave, taphTrend, telvet]

Return type str

anesplot.record_main.**plot_trenddata**(datadf, header, param_dico)

clinical main plots of a trend recordings

parameters df : pdDataframe

recorded data (MonitorTrend.data)

header [dict] recording parameters (MonitorTrend.header)

param_dico [dict] plotting parameters (MonitorTrend.param)

Returns **afig_dico**

Return type dict of name:fig

class anesplot.record_main.**MonitorTrend**(*filename=None, load=True*)

Bases: anesplot.record_main._SlowWave

monitor trends recordings:

input = filename : path to file load = boolean to load data (default is True)

file [str] short name

filename [str] long name

header [dict] record parameters

param [dict] parameters

clean_trend [external] clean the data

show_graphs [external] plot clinical main plots

class anesplot.record_main.**TaphTrend**(*filename=None*)

Bases: anesplot.record_main._SlowWave

taphonius trends recordings

input ... FILLME

attributes ... FILLME

extract_taph_actions(*data*)

extract Taph actions

Parameters **data** (*pandas dataframe*) – record df form taphonius recording)

Returns

Return type actiondf pandas dataframe

class anesplot.record_main.**TelevetWave**(*filename=None*)

Bases: anesplot.record_main._FastWave

class to organise teleVet recordings transformed to csv files. input:

filename : str (fullpath, default:None)

class anesplot.record_main.**MonitorWave**(*filename=None, load=True*)

Bases: anesplot.record_main._FastWave

class to organise monitorWave recordings. input : filename = path to file load = boolean to load data (default is True)

attibutes ... FILLME

methods ... FILLME

anesplot.record_main.**main**(*file_name=None*)

main script called from command line call : “python anesthPlot/anesplot/__main__.py” args : optional filename (fullname)

return: set of plots for either monitorTrend, monitorWave oe televet recording

MODULES

3.1 anesplot package

3.1.1 Subpackages

`anesplot.config` package

Submodules

`anesplot.config.build_recordRc` module

`anesplot.config.load_recordRc` module

Module contents

`anesplot.loadrec` package

Submodules

`anesplot.loadrec.explore` module

Created on Thu Mar 12 16:52:13 2020

@author: cdesbois

`anesplot.loadrec.explore.gui_choosefile(paths=None)`
select a file via a dialog and return the file name.

`anesplot.loadrec.loadmonitor_trendrecord` module

Created on Wed Jul 24 13:43:26 2019 @author: cdesbois

load a monitor trend recording:

- choose a file
- load the header to a dictionary
- load the date into a pandas dataframe

`anesplot.loadrec.loadmonitor_trendrecord.choosefile_gui(dirname=None)`

Select a file via a dialog and return the (full) filename.

Parameters `dir_path` (*str*) – location to place the gui ('generally paths['data']) else home

Returns `fname[0]` – filename

Return type `str`

`anesplot.loadrec.loadmonitor_trendrecord.loadmonitor_trendheader(filename)`

load the file header.

Parameters `filename` (*str*) – full name of the file

Returns `header`

Return type `dict`

`anesplot.loadrec.loadmonitor_trendrecord.loadmonitor_trenddata(filename, headerdico)`

load the monitor trend data

Parameters

- `filename` (*str*) – fullname
- `headerdico` (*dict*) – fileheader

Returns `df` = trends data

Return type `pandas.DataFrame`

anesplot.loadrec.loadmonitor_waverecord module

Created on Wed Jul 24 14:56:58 2019 @author: cdesbois

load a monitor wave recording:

- choose a file
- load the header to a pandas dataframe
- load the data into a pandas dataframe

`anesplot.loadrec.loadmonitor_waverecord.choosefile_gui(dirname=None)`

Select a file via a dialog and return the (full) filename.

Parameters `dir_path` (*str*) – location to place the gui ('generally paths['data']) else home

Returns `fname[0]` – filename

Return type `str`

`anesplot.loadrec.loadmonitor_waverecord.loadmonitor_waveheader(filename=None) → dict`

load the wave file header.

Parameters `filename` (*str*) – full name of the file

Returns `header`

Return type `dictionary`

`anesplot.loadrec.loadmonitor_waverecord.loadmonitor_wavedata(filename=None)`
load the monitor wave csvDataFile.

Parameters `filename` (*str*) – full name of the file

Returns `df` = trends data

Return type `pandas.DataFrame`

anesplot.loadrec.loadtaph_trendrecord module

Created on Wed Jul 24 15:30:07 2019 @author: cdesbois

load a taphonius data recording:

- choose a file
- load the patient datafile to a dictionary
- load the physiological date into a pandas dataframe

nb = 4 files per recording :

- .pdf -> anesthesia record 'manual style'
- .xml -> taphonius technical record -> to be extracted
- Patient.csv -> patient id and specifications
- SD...csv -> anesthesia record

`anesplot.loadrec.loadtaph_trendrecord.build_taph_decodedate_dico(pathdict=None)`
list all the taph recordings and the paths to the record: input:
paths: dictionary containing { 'taph': pathToTheData }

output: dictionary: {date : filename}

`anesplot.loadrec.loadtaph_trendrecord.extract_record_day(monitor_file_name)`
extract the date as 'YYYY_MM_DD' from a monitor_filename input:
monitor file name (shortname)

output: day : YYYY_MM_DD str

`anesplot.loadrec.loadtaph_trendrecord.choose_taph_record(monitorname=None)`
select the taph recording: input:
taphdico : {date:path} builded from build_taph_decodedate_dico() 'year = integer to place the pointer in pull down menu date = to be implemented (as year but to extract from monitor filename)

output: filename (str) full path

`anesplot.loadrec.loadtaph_trendrecord.loadtaph_trenddata(filename)`
load the taphoniusData trends data.

Parameters `filename` (*str*) – fullname

Returns `df` = trends data

Return type `pandas.DataFrame`

`anesplot.loadrec.loadtaph_trendrecord.loadtaph_patientfile(filename)`

load the taphonius patient.csv file input:

filename [(str) the full filename] (the headername will be reconstructed inside the function)

output: descr = dict of patient_data

anesplot.loadrec.loadtelevet module

Created on Wed Jul 31 16:22:06 2019 @author: cdesbois

load televet exported (csv) data: to be developped

`anesplot.loadrec.loadtelevet.choosefile_gui(dirpath=None)`

select a file using a dialog.

Parameters **dir_path** (str) – optional location of the data (ex : paths[‘data’], default : ‘home’)

Returns filename (full path)

Return type str

`anesplot.loadrec.loadtelevet.loadtelevet(fname=None, all_traces=False)`

load the televetCsvExportedFile.

Parameters

- **file** (str) – name of the file
- **all_traces** (bool) – load all the derivations

Returns df = recorded traces

Return type pandas.DataFrame

Module contents

anesplot.plot package

Submodules

anesplot.plot.trend_plot module

Created on Tue Apr 19 09:08:56 2016 @author: cdesbois

collection of functions to plot the trend data

`anesplot.plot.trend_plot.remove_outliers(df, key, limits=None)`

remove outliers input:

df : pandas.DataFrame key : a column label limits : dictionary of key:(limLow, limHigh)

output: pandas.series without the outliers

`anesplot.plot.trend_plot.color_axis(ax, spine='bottom', color='r')`
change the color of the label & tick & spine.

Parameters

- **ax** (*matplotlib.pyplot.axis*) – the axis
- **spine** (*str*) – optional location in ['bottom', 'left', 'top', 'right']
- **colors** (*str*) – optional color

`anesplot.plot.trend_plot.append_loc_to_fig(ax, dt_list, label='g')`
append vertical lines to indicate a location 'eg: arterial blood gas'

Parameters

- **ax** (*matplotlib.pyplot.axis*) – the axis
- **dt_list** (*[datetime]*) – list of datetime values
- **label** (*str*) – a key to add to the label (default is 'g')

Returns **res** a dictionary containing the locations

Return type dict

`anesplot.plot.trend_plot.save_graph(path, ext='png', close=True, verbose=True)`
Save a figure from pyplot. :param path: The path (and filename, without the extension) to save the figure to.

Parameters

- **ext** (*string (default='png')*) – The file extension. This must be supported by the active matplotlib backend (see matplotlib.backends module). Most backends support 'png', 'pdf', 'ps', 'eps', and 'svg'.
- **close** (*boolean (default=True)*) – Whether to close the figure after saving. If you want to save the figure multiple times (e.g., to multiple formats), you should NOT close it in between saves or you will have to re-plot it.
- **verbose** (*boolean (default=True)*) – Whether to print information about when and where the image has been saved.

`anesplot.plot.trend_plot.plot_header(descr, param=None)`
plot the header of the file.

Parameters

- **descr** (*dict*) – header of the recording
- **param** (*dict*) – dictionary of parameters

Returns **fig** plot of the header

Return type pyplot.figure

`anesplot.plot.trend_plot.hist_cardio(data, param=None)`
mean arterial pressure histogramme using matplotlib.

Parameters

- **data** (*pandas.DataFrame*) – the recorded trends data (keys used : 'ip1m' and 'hr'),
- **param** (*dict*) – parameters (save=boolean, 'path': path to directory)

Returns **fig** matplotlib.pyplot.figure

`anesplot.plot.trend_plot.plot_one_over_time(x, y, colour)`
 plot y over x using colour

`anesplot.plot.trend_plot.hist_co2_iso(data, param=None)`
 CO2 and iso histogramme (NB CO2 should have been converted from % to mmHg)

Parameters

- **data** (*pandas.DataFrame*) – the trends recorded data
- **param** (*dict*) – dictionary of parameters

Returns fig pyplot.figure

`anesplot.plot.trend_plot.cardiovasc(data, param=None)`
 cardiovascular plot

Parameters

- **data** (*pandas.DataFrame*) – the recorded trends data keys used :['ip1s', 'ip1m', 'ip1d', 'hr']
- **param** (*dict*) – dict(save: boolean, path['save'], xmin, xmax, unit, dtime = boolean for time display in HH:MM format)

Returns fig= pyplot.figure

`anesplot.plot.trend_plot.cardiovasc_p1p2(data, param=None)`
 cardiovascular plot with central venous pressure (p2)

Parameters

- **data** (*pandas.DataFrame*) – the trends recorded data keys used :['ip1s', 'ip1m', 'ip1d', 'hr', 'ip2s', 'ip2m', 'ip2d']
- **param** (*dict*) – dict(save: boolean, path['save'], xmin, xmax, unit, dtime = boolean for time display in HH:MM format)

Returns fig= pyplot.figure

`anesplot.plot.trend_plot.co2iso(data, param=None)`
 anesth plot (CO2/iso)

Parameters

- **data** (*pandas.DataFrame*) – the recorded data keys used :['ip1s', 'ip1m', 'ip1d', 'hr']
- **param** (*dictionary*) – dict(save: boolean, path['save'], xmin, xmax, unit, dtime = boolean for time display in HH:MM format)

:returns fig= pyplot.figure

`anesplot.plot.trend_plot.func(ax, x, y1, y2, color='tab:blue', x0=38)`

`anesplot.plot.trend_plot.co2o2(data, param)`
 respiratory plot (CO2 and Iso)

Parameters

- **data** (*pandas.DataFrame*) – recorded trends data keys used :['ip1s', 'ip1m', 'ip1d', 'hr']
- **param** (*dict*) – dict(save: boolean, path['save'], xmin, xmax, unit, dtime = boolean for time display in HH:MM format)

Returns fig= pyplot.figure

`anesplot.plot.trend_plot.ventil(data, param)`
 plot ventilation parameters (.tvInsp, .pPeak, .pPlat, .peep, .minVexp, .co2RR, .co2exp)

Parameters

- **data** (*pandas.DataFrame*) – recorded data, keys used :['ip1s', 'ip1m', 'ip1d', 'hr']
- **param** (*dict*) – dict(save: boolean, path['save'], xmin, xmax, unit, dtime = boolean for time display in HH:MM format)

Returns fig= pyplot.figure

`anesplot.plot.trend_plot.recruit(data, param)`
 display a recruit manoeuver (.pPeak, .pPlat, .peep, .tvInsp)

Parameters

- **data** (*pandas.DataFrame*) – recorded data keys used :['ip1s', 'ip1m', 'ip1d', 'hr']
- **param** (*dict*) – dict(save: boolean, path['save'], xmin, xmax, unit, dtime = boolean for time display in HH:MM format)

:returns fig= pyplot.figure

`anesplot.plot.trend_plot.ventil_cardio(data, param)`
 build ventilation and cardiovascular plot

Parameters

- **data** (*pandas.DataFrame*) – teh recorded trends data keys used :['ip1s', 'ip1m', 'ip1d', 'hr']
- **param** (*dict*) – dict(save: boolean, path['save'], xmin, xmax, unit, dtime = boolean for time display in HH:MM format)

Returns fig= pyplot.figure

`anesplot.plot.trend_plot.save_distri(data, path)`
 save as 'O_.' the 4 distributions graphs for cardiovasc annd respi

`anesplot.plot.trend_plot.fig_memo(path, fig_name)`
 append latex citation commands in a txt file inside the fig folder create the file iif it doesn't exist

anesplot.plot.wave_plot module

Created on Tue Apr 19 09:08:56 2016

@author: cdesbois

`anesplot.plot.wave_plot.color_axis(ax, spine='bottom', color='r')`
 change the color of the label & tick & spine.

Parameters

- **ax** (*matplotlib.pyplot.axis*) – the axis
- **spine** (*str*) – optional location in ['bottom', 'left', 'top', 'right']
- **colors** (*str*) – optional color

`anesplot.plot.wave_plot.plot_wave(data, keys, param)`
 plot the waves recorded (from as5)

Parameters

- **data** (*pandas.DataFrame*) – the recorded trends data
- **keys** (*list*) – one or two in ['wekg','ECG','wco2','wawp','wflow','wap']
- **{mini}** (*dict*) – limits in point value (index), maxi: limits in point value (index)}

Returns fig plt.figure the plot

Returns lines plt.line2D the line to animate

(Nb plot data/index, but the xscale is indicated as sec)

`anesplot.plot.wave_plot.get_roi(waves)`

use the drawn figure to extract the relevant data in order to build an animation

Parameters waves (*MonitorWave object*) – a wave recording

Returns a dictionary containing ylims, xlims(point, dtime and sec),

traces used to build the plot, the fig object :rtype: dictionary

`anesplot.plot.wave_plot.create_video(waves, speed=1, save=False, savename='example', savedir='~')`

create a video from a figure input:

waves : waves object speed : integer, speed of the display save : boolean (default=False) savename : str (default='example') savedir : str (path, default='~')

Returns .mp4 file .png file

Module contents

Created on Tue Apr 19 09:08:56 2016

functions to plot the trend data

@author: cdesbois

anesplot.treatrec package

Submodules

anesplot.treatrec.clean_data module

Created on Wed Jul 31 16:05:29 2019

@author: cdesbois

`anesplot.treatrec.clean_data.clean_trenddata(df)`

remove artifacts in the recorded trends

anesplot.treatrec.ekg_to_hr module

Created on Wed Feb 12 16:52:00 2020 @author: cdesbois

function used to treat an EKG signal and extract the heart rate typically (copy, paste and execute line by line)

0. after

```
:: import pandas as pd
import anesplot.record_main as rec from anesplot.treatrec import ekg_to_hr as tohr
```

1. load the data in a pandas dataframe:

(through classes rec.MonitorTrend & rec.MonitorWave)

```
trendname = '' # fullname
or
trendname = rec.choosefile_gui()
```

```
wavename = rec.trendname_to_wavename(trendname)
-
# load the data
trends = rec.MonitorTrend(trendname)
waves = rec.MonitorWave(wavename)
-
# format the name
name = trends.header['Patient Name'].title().replace(' ', '')
name = name[0].lower() + name[1:]
```

2. treat the ekg wave:

- get parameters
- build a dataframe to work with (waves)
- low pass filtering
- build the beat locations (beat based dataframe):

```
params = waves.param
ekg_df = pd.DataFrame(waves.data.wekg)
ekg_df['wekg_lowpass'] = rec.wf.fix_baseline_wander(ekg_df.wekg,
                                                    waves.param['fs'])
beat_df = tohr.detect_beats(ekg_df.wekg_lowpass, mult=1)
```

3. perform the manual adjustments required:

- based on a graphical display of beat locations, an rr values
- build a container for the manual corrections:

```
figure = tohr.plot_beats(ekg_df.wekg_lowpass, beat_df)
to_change_df = pd.DataFrame(columns=beat_df.columns.insert(0, 'action'))
```

- remove or add peaks : zoom on the figure to observe only one peak, then:

```
to_change_df = tohr.remove_beat(beat_df, ekg_df, to_change_df, figure)
or
to_change_df = tohr.append_beat(beat_df, ekg_df, to_change_df, figure,
                               yscale=1)
```

- combine to update the beat_df with the manual changes:

```
beat_df = tohr.update_beat_df(beat_df, to_change_df,
                              path_to_file="", from_file=False)
```

- save the peaks locations:

```
tohr.save_beats(beat_df, to_change_df, savename='', savepath=None)
(# or reload
beat_df = pd.read_hdf('beatDf.hdf', key='beatDf') )
```

4. go from points values to continuous time:

```
beat_df = tohr.compute_rr(beat_df)
ahr_df = tohr.interpolate_rr(beat_df)
tohr.plot_rr(ahr_df, params)
```

5. append intantaneous heart rate to the initial data:

```
ekg_df = tohr.append_rr_and_ihr_to_wave(ekg_df, ahr_df)
waves.data = tohr.append_rr_and_ihr_to_wave(waves.data, ahr_df)
trends.data = tohr.append_ihr_to_trend(trends.data, waves.data, ekg_df)
```

6. save:

```
tohr.save_trends_data(trends.data, savename=name, savepath='data')
tohr.save_waves_data(waves.data, savename=name, savepath='data')
```

```
anesplot.treatrec.ekg_to_hr.detect_beats(ser, fs=300, species='horse', mult=1)
```

detect the peak locations

Parameters

- **ser** (*pandas.series*) – the data
- **fs** (*integer*) – sampling frequency
- **species** (*string*) – in [horse]
- **mult** (*float*) – correction / 1 for qRs amplitude

Returns df=pandas.DataFrame

anesplot.treatrec.ekg_to_hr.**plot_beats**(*ecg, beats*)
plot ecg waveform + beat location

anesplot.treatrec.ekg_to_hr.**append_beat**(*beatdf, ekgdf, tochange_df, fig, lim=None, yscale=1*)
locate the beat in the figure, append to a dataframe['toAppend']

Parameters

- **beatdf** (*pandas.DataFrame*) – contains the point based location (pLocs)
- **ekgdf** (*pandas dataframe*) – contains the wave recording ((wekg_lowpass)
- **tochange_df** (*pandas.DataFrame*) – to store the beats toAppend or toRemove
- **fig** (*pyplot.Figure*) – figure to find time limits
- **lim** (*integer*) – ptBasedLim optional to give it manually
- **yscale** (*float*) – amplitude mutliplication factor for detection (default=1)

Returns tochange_df: incremented changedf (pt location)

Return type pandasDataFrame

methods :

locate the beat in the figure, append to a dataframe['toAppend'] 0.: if not present : build a dataframe:

```
>>> to_change_df = pd.DataFrame(columns=['toAppend', 'toRemove'])
```

1.: locate the extra beat in the figure (cf plot_beats()) and zoom to observe only a negative peak

2.: call the function:

```
>>> to_change_df = remove_beat(beatdf, ekgdf, tochange_df, fig)
-> the beat parameters will be added the dataframe
```

.in the end of the manual check, update the beat_df

- **first** : save beat_df and to_change_df
- **second** [run:]

```
>>> beat_df = update_beat_df()
```

anesplot.treatrec.ekg_to_hr.**remove_beat**(*beatdf, ekgdf, tochange_df, fig, lim=None*)
locate the beat in the figure, append to a dataframe['toRemove']

0.: if not present build a dataframe:

```
>>> to_change_df = pd.DataFrame(columns=['toAppend', 'toRemove'])
```

1.: locate the extra beat in the figure (cf plot_beats()) and zoom to observe only a negative peak

2.: call the function:::

```
>>> to_change_df = remove_beat(beatdf, ekgdf, tochange_df, fig)
-> the beat parameters will be added the DataFrame
```

.(in the end of the manual check, update the beat_df

- first : save beat_df and to_change_df
- **second** [run]

```
>>> beat_df = update_beat_df()
```

`anesplot.treatrec.ekg_to_hr.save_beats(beatdf, tochangedf, savename="", dirpath=None)`
 save the beats locations as csv and hdf5 file

Parameters

- **beatde** (*pd.dataframes*) –
- **tochangedf** (*pandas.dataframe*) –
- **savename** (*filename*) –
- **dirpath** (*path to save in*) –
- **output** –
- -----
- **file** (*hdf*) –
- **key='beatDf'** –

`anesplot.treatrec.ekg_to_hr.update_beat_df(beatdf, tochangedf, path_to_file="", from_file=False)`
 implement in the beat location the manual corrections fromFile = True force the disk loading of the dataframes

`anesplot.treatrec.ekg_to_hr.compute_rr(beatdf, fs=None)`
 compute rr intervals (from pt to time)

Parameters

- **beatdf** (*pd.DataFrame*) – with 'pLoc'
- **fs** (*integer*) – sampling frequency

Returns with: 'rr' = rr duration 'rrDiff' = rrVariation 'rrSqDiff' = rrVariation^2

Return type *pd.DataFrame*

`anesplot.treatrec.ekg_to_hr.interpolate_rr(beatdf, kind=None)`
 interpolate the beat_df (pt -> time values)

Parameters

- **beatDf** (*pd.DataFrame*) –
- **kind** (*str*) – 'linear' or 'cubic'(default)

Returns 'espts' = evenly spaced points 'rrInterpol' = interpolated rr

Return type *pd.DataFrame* with evenly spaced data

`anesplot.treatrec.ekg_to_hr.plot_rr(ahr_df, param, HR=False)`
 plot RR vs pt values + rrSqDiff

Parameters

- **pdDataFrame** (*hr_df* =) –

- **params** – dict containing 'fs' as key

`anesplot.treatrec.ekg_to_hr.append_rr_and_ihr_to_wave(wave, ahrdf)`
 append rr and ihr to the waves based on pt value (ie index)

`anesplot.treatrec.ekg_to_hr.plot_agreement(trenddf)`
 plot ip1HR & ihr to check agreement

`anesplot.treatrec.ekg_to_hr.append_ihr_to_trend(trenddf, wavedf, ekgdf)`
 append 'ihr' (instantaneous heart rate) to the trends

`anesplot.treatrec.ekg_to_hr.save_trends_data(trenddf, savename='', dirpath='data')`

save the trends data to a csv and hd5 file, including an ihr column

trenddf : pd.dataframes savename : str dirpath : str

path to save in (default= current working directory)

hdf file, key='trends_data'

`anesplot.treatrec.ekg_to_hr.save_waves_data(wavedf, savename='', dirpath='data')`
 save the trends data to a hd5 file, including an ihr column

Parameters

- **trenddf** (*pd.dataframes*) –
- **savename** (*str*) – dirpath : path to save in (default='data')
- **output** –
- ----- –
- **hdf_file** –
- **key='waves_data'** –

anesplot.treatrec.extract_hypotension module

Spyder Editor

This is a temporary script file.

`anesplot.treatrec.extract_hypotension.extract_hypotension(atrend, pamin=70)`
 return a dataframe with the beginning and ending phses of hypotension

Parameters

- **atrend** (*MonitorTrend object*) –
- **pamin** (*float= threshold de define hypotension on mean arterial pressure*) –
- **70** (*((default is)*) –

Returns **durdf** – transitions (up and down, in seconds from beginning) and duration in the hypotension state (in seconds)

Return type pandas DataFrame containing

`anesplot.treatrec.extract_hypotension.plot_hypotension(atrend, durdf, durmin=15, pamin=70)`
 plot the hupotentions phases

Parameters

- **atrend** (*TYPE*) – DESCRIPTION.
- **durdf** (*TYPE*) – DESCRIPTION.
- **durmin** (*TYPE*, *optional*) – DESCRIPTION. The default is 15.

Returns *fig* – DESCRIPTION.

Return type *TYPE*

`anesplot.treatrec.extract_hypotension.scatter_length_meanhypo(atrend, durdf)`
 draw a scatter plot (hypotensive arterial value vs duration of hypotension) :param trends: :type trends: MonitorTrend :param durdf: :type durdf: pandas dataframe containing the value and duration

Returns *fig*

Return type matplotlib.pyplot figure

`anesplot.treatrec.extract_hypotension.plot_all_dir_hypo(dirname=None, scatter=False)`
 walk through the folder and plot the values

anesplot.treatrec.hr_to_hrv module

`anesplot.treatrec.hr_to_hrv.build_hrv_limits(spec='horse')`
 return a dico containing HRV limits (VLF, LF, HF) input : spec in ['horse', 'man']

anesplot.treatrec.wave_func module

Created on Fri Dec 8 12:46:41 2017

@author: cdesbois

`anesplot.treatrec.wave_func.fix_baseline_wander(data, fs=500)`
 BaselineWanderRemovalMedian.m from ecg-kit. Given a list of amplitude values (data) and sample rate (sr), it applies two median filters to data to compute the baseline. The returned result is the original data minus this computed baseline.

`anesplot.treatrec.wave_func.rol_mean(ser, win_lengh=1, fs=500)`
 returns a rolling mean of a RR serie

Parameters

- **pd.Series** (*ser*) –
- **win_lengh** (*integer*) – window lenght for averaging (in sec),
- **fs** (*int*) – sampling frequency

`anesplot.treatrec.wave_func.return_points(df, fig)`
 return a tuple containing the point values of ROI

Parameters

- **df** (*anesthesia record dataframe*) –
- **fig** (*pyplot.figure*) –

Returns ROI

Return type dict

`anesplot.treatrec.wave_func.restrict_time_area(df1, mini=None, maxi=None)`
 return a new dataframe with reindexation

Parameters

- **df1** (*pandas.DataFrame*) –
- **mini** (*integer*) – miniPointValue
- **maxi** (*integer*) – maxiPointValue

Returns

Return type `pandas.DataFrame`

Module contents**3.1.2 Submodules****3.1.3 Module contents**

anesthPlot is a package to plot/use clinical anesthesia records for teaching

three way to use it:

1. **run directly anesplot from a terminal** -> `PYTHONPATH=<pathToAnesthPlot> python -m anesplot ->`
generate a quick plotting of most interesting parts
2. **from an ipython terminal** -> `import anesthPlot.anesplot.recordmain as rec -> trends =`
`rec.MonitorTrend() -> waves = rec.MonitorWave() -> ...` and use the objects trends and waves
3. import the module in a python environment (see below)

(the presets are actually designed

- for use with equine anesthesia
- to load data from a Monitor generated datex AS3/5 monitoring machine)

typical use when importing the module to build a clinical case

```
import os
import sys
```

```
import numpy as np
import pandas as pd
```

```
import anesthPlot.record_main as rec
sys.path.append(os.path.expanduser('~/.pg/Utils'))
from Utils import saveGraph
import bloodGases.bgmain_manual as bgman
```

```
paths = rec.paths
paths['save'] = os.path.expanduser('~/.toPlay/temp/')
os.chdir(paths['save'])
```

```
## globals
def save_plot(name):
```

```
    filename = os.path.join(paths['save'], 'fig', name)
    saveGraph(filename, ext='png', close=False, verbose=True)
```

def explore_hdf(filename):

```
    try:
        hdf = pd.HDFStore(filename)
        keys = [key.replace('/', '') for key in hdf.keys()]
        print(' found h5_file { } that contains { } '.format(filename, keys))
```

```
        hdf.close()
```

```
    except:
        print('{ } is not an h5 file'.format(filename))
```

```
saveName = os.path.join(paths['save'], 'data', 'aname.h5')
explore_hdf(saveName)
## load and work trendName = rec.choosefile_gui(paths['data']) WaveName = rec.choosefile_gui(paths['data'])
# build objects with headers trends = rec.MonitorTrend(trendName, load=True) waves = rec.MonitorWave(waveName,
load=True
# or append data (pretreated ones) #trends.data = pd.read_hdf(saveName, 'trend_df') #waves.data =
pd.read_hdf(saveName, 'wave_df')
#remove filenames del waveName, trendName
# now you are ready to work with loaded trends and waves
```


INDICES AND TABLES

- `genindex`
- `modindex`
- `search`

PYTHON MODULE INDEX

a

- `anesplot`, 19
- `anesplot.config`, 6
- `anesplot.loadrec`, 9
 - `anesplot.loadrec.explore`, 6
 - `anesplot.loadrec.loadmonitor_trendrecord`, 6
 - `anesplot.loadrec.loadmonitor_waverecord`, 7
 - `anesplot.loadrec.loadtaph_trendrecord`, 7
 - `anesplot.loadrec.loadtelevet`, 8
- `anesplot.plot`, 13
 - `anesplot.plot.trend_plot`, 9
 - `anesplot.plot.wave_plot`, 12
- `anesplot.record_main`, 3
- `anesplot.treatrec`, 19
 - `anesplot.treatrec.clean_data`, 13
 - `anesplot.treatrec.ekg_to_hr`, 13
 - `anesplot.treatrec.extract_hypotension`, 18
 - `anesplot.treatrec.hr_to_hrv`, 18
 - `anesplot.treatrec.wave_func`, 19

A

[adapt_with_syspath\(\)](#) (in module *anesplot.config.load_recordrc*), 6
[anesplot](#)
 module, 19
[anesplot.config](#)
 module, 6
[anesplot.config.build_recordrc](#)
 module, 5
[anesplot.config.load_recordrc](#)
 module, 5
[anesplot.loadrec](#)
 module, 9
[anesplot.loadrec.explore](#)
 module, 6
[anesplot.loadrec.loadmonitor_trendrecord](#)
 module, 6
[anesplot.loadrec.loadmonitor_waverecord](#)
 module, 7
[anesplot.loadrec.loadtaph_trendrecord](#)
 module, 7
[anesplot.loadrec.loadtelevet](#)
 module, 8
[anesplot.plot](#)
 module, 13
[anesplot.plot.trend_plot](#)
 module, 9
[anesplot.plot.wave_plot](#)
 module, 12
[anesplot.record_main](#)
 module, 3
[anesplot.treatrec](#)
 module, 19
[anesplot.treatrec.clean_data](#)
 module, 13
[anesplot.treatrec.ekg_to_hr](#)
 module, 13
[anesplot.treatrec.extract_hypotension](#)
 module, 18
[anesplot.treatrec.hr_to_hrv](#)
 module, 18
[anesplot.treatrec.wave_func](#)

module, 19

[append_beat\(\)](#) (in module *anesplot.treatrec.ekg_to_hr*), 15
[append_ihr_to_trend\(\)](#) (in module *anesplot.treatrec.ekg_to_hr*), 17
[append_loc_to_fig\(\)](#) (in module *anesplot.plot.trend_plot*), 9
[append_rr_and_ihr_to_wave\(\)](#) (in module *anesplot.treatrec.ekg_to_hr*), 17

B

[build_hrv_limits\(\)](#) (in module *anesplot.treatrec.hr_to_hrv*), 18
[build_paths\(\)](#) (in module *anesplot.config.load_recordrc*), 5
[build_taph_decodedate_dico\(\)](#) (in module *anesplot.loadrec.loadtaph_trendrecord*), 8

C

[cardiovasc\(\)](#) (in module *anesplot.plot.trend_plot*), 10
[cardiovasc_p1p2\(\)](#) (in module *anesplot.plot.trend_plot*), 10
[choose_taph_record\(\)](#) (in module *anesplot.loadrec.loadtaph_trendrecord*), 8
[choosefile_gui\(\)](#) (in module *anesplot.loadrec.loadmonitor_trendrecord*), 6
[choosefile_gui\(\)](#) (in module *anesplot.loadrec.loadmonitor_waverecord*), 7
[choosefile_gui\(\)](#) (in module *anesplot.loadrec.loadtelevet*), 8
[choosefile_gui\(\)](#) (in module *anesplot.record_main*), 3
[clean_trenddata\(\)](#) (in module *anesplot.treatrec.clean_data*), 13
[co2iso\(\)](#) (in module *anesplot.plot.trend_plot*), 11
[co2o2\(\)](#) (in module *anesplot.plot.trend_plot*), 11
[color_axis\(\)](#) (in module *anesplot.plot.trend_plot*), 9
[color_axis\(\)](#) (in module *anesplot.plot.wave_plot*), 12
[compute_rr\(\)](#) (in module *anesplot.treatrec.ekg_to_hr*), 16
[create_video\(\)](#) (in module *anesplot.plot.wave_plot*), 12

D

`detect_beats()` (in module *anesplot.treatrec.ekg_to_hr*), 15

E

`extract_hypotension()` (in module *anesplot.treatrec.extract_hypotension*), 18

`extract_record_day()` (in module *anesplot.loadrec.loadtaph_trendrecord*), 8

`extract_taph_actions()` (*anesplot.record_main.TaphTrend* method), 4

F

`fig_memo()` (in module *anesplot.plot.trend_plot*), 12

`filedialog()` (in module *anesplot.config.build_recordrc*), 5

`fix_baseline_wander()` (in module *anesplot.treatrec.wave_func*), 19

`func()` (in module *anesplot.plot.trend_plot*), 11

G

`get_roi()` (in module *anesplot.plot.wave_plot*), 12

`gui_choosefile()` (in module *anesplot.loadrec.explore*), 6

H

`hist_cardio()` (in module *anesplot.plot.trend_plot*), 10

`hist_co2_iso()` (in module *anesplot.plot.trend_plot*), 10

I

`interpolate_rr()` (in module *anesplot.treatrec.ekg_to_hr*), 17

L

`loadmonitor_trenddata()` (in module *anesplot.loadrec.loadmonitor_trendrecord*), 6

`loadmonitor_trendheader()` (in module *anesplot.loadrec.loadmonitor_trendrecord*), 6

`loadmonitor_wavedata()` (in module *anesplot.loadrec.loadmonitor_waverecord*), 7

`loadmonitor_waveheader()` (in module *anesplot.loadrec.loadmonitor_waverecord*), 7

`loadtaph_patientfile()` (in module *anesplot.loadrec.loadtaph_trendrecord*), 8

`loadtaph_trenddata()` (in module *anesplot.loadrec.loadtaph_trendrecord*), 8

`loadtelevet()` (in module *anesplot.loadrec.loadtelevet*), 8

M

`main()` (in module *anesplot.config.build_recordrc*), 5

`main()` (in module *anesplot.record_main*), 4

module

anesplot, 19

anesplot.config, 6

anesplot.config.build_recordrc, 5

anesplot.config.load_recordrc, 5

anesplot.loadrec, 9

anesplot.loadrec.explore, 6

anesplot.loadrec.loadmonitor_trendrecord, 6

anesplot.loadrec.loadmonitor_waverecord, 7

anesplot.loadrec.loadtaph_trendrecord, 7

anesplot.loadrec.loadtelevet, 8

anesplot.plot, 13

anesplot.plot.trend_plot, 9

anesplot.plot.wave_plot, 12

anesplot.record_main, 3

anesplot.treatrec, 19

anesplot.treatrec.clean_data, 13

anesplot.treatrec.ekg_to_hr, 13

anesplot.treatrec.extract_hypotension, 18

anesplot.treatrec.hr_to_hrv, 18

anesplot.treatrec.wave_func, 19

MonitorTrend (class in *anesplot.record_main*), 4

MonitorWave (class in *anesplot.record_main*), 4

P

`plot_agreement()` (in module *anesplot.treatrec.ekg_to_hr*), 17

`plot_all_dir_hypo()` (in module *anesplot.treatrec.extract_hypotension*), 18

`plot_beats()` (in module *anesplot.treatrec.ekg_to_hr*), 15

`plot_header()` (in module *anesplot.plot.trend_plot*), 10

`plot_hypotension()` (in module *anesplot.treatrec.extract_hypotension*), 18

`plot_one_over_time()` (in module *anesplot.plot.trend_plot*), 10

`plot_rr()` (in module *anesplot.treatrec.ekg_to_hr*), 17

`plot_trenddata()` (in module *anesplot.record_main*), 3

`plot_wave()` (in module *anesplot.plot.wave_plot*), 12

R

`read_config()` (in module *anesplot.config.build_recordrc*), 5

`recrut()` (in module *anesplot.plot.trend_plot*), 11

`remove_beat()` (in module *anesplot.treatrec.ekg_to_hr*), 16

`remove_outliers()` (in module *anesplot.plot.trend_plot*), 9

`restrict_time_area()` (in module *anesplot.treatrec.wave_func*), 19

`return_points()` (in module *anesplot.treatrec.wave_func*), 19

`rol_mean()` (in module *anesplot.treatrec.wave_func*), 19

S

`save_beats()` (in module *anesplot.treatrec.ekg_to_hr*), 16

`save_distri()` (in module *anesplot.plot.trend_plot*), 12

`save_graph()` (in module *anesplot.plot.trend_plot*), 9

`save_trends_data()` (in module *anesplot.treatrec.ekg_to_hr*), 17

`save_waves_data()` (in module *anesplot.treatrec.ekg_to_hr*), 17

`scatter_length_meanhypo()` (in module *anesplot.treatrec.extract_hypotension*), 18

`select_type()` (in module *anesplot.record_main*), 3

`select_wave()` (in module *anesplot.record_main*), 3

T

`TaphTrend` (class in *anesplot.record_main*), 4

`TelevetWave` (class in *anesplot.record_main*), 4

`trendname_to_wavename()` (in module *anesplot.record_main*), 3

U

`update_beat_df()` (in module *anesplot.treatrec.ekg_to_hr*), 16

V

`ventil()` (in module *anesplot.plot.trend_plot*), 11

`ventil_cardio()` (in module *anesplot.plot.trend_plot*), 11

W

`write_configfile()` (in module *anesplot.config.build_recordrc*), 5