
anesthPlot

Release beta

Christophe Desbois

Sep 11, 2021

MAIN_{*s*}CRIP*T* :

1	Welcome to anesthPlot's documentation!	1
2	main script	3
3	modules	7
4	Indices and tables	21
	Python Module Index	23

WELCOME TO ANESTHPLOT'S DOCUMENTATION!

anesthPlot is a python package developped to extract, manipulate and plots anesthesia data recorded from the Monitor Software to be used mostly in a teaching environment.

Warning: This project is:

- a work in process
- the processes are mainly focused on horses anesthesia
- in our environment the data recorded came from an as3 or as5 anesthesia machine

1.1 Features

- **load** recordings from a trend or a wave recordings
- build a **standard debriefing** (trends) **plot series** (script usage)
 - global histograms (cardiovascular and anesthesia summary)
 - cardiovascular trends time based plots
 - respiratory trends time based plots
 - anesthesia trends time based plots
- build a **plot for wave** recording (one or two waves (script usage)
- can be used as a **python package**
 - **usage :**

```
import anesthplot.record_main as rec
trendname = 'a_full_path_to_csv_file'
trends = rec.MonitorTrend(trendname)
wavename = rec.trendname_to_wavename(trendname)
waves = rec.MonitorWave(trends)
trends.show_graphs()
```


MAIN SCRIPT

2.1 anesplot.record_main module

main script/module to load and display an anesthesia record

can be runned as a script:: python record_main.py

or imported as a package:: import anesplot.record_main as rec

`anesplot.record_main.choosefile_gui(dir_path=None)`

Select a file via a dialog and return the (full) filename.

Parameters `dir_path` (*str*) – location to place the gui ('generally paths['data']) else home

Returns `fname[0]` – filename

Return type *str*

`anesplot.record_main.trendname_to_wavename(name)`

just compute the supposed name

`anesplot.record_main.select_type(question=None, items=None, num=0)`

select the recording type:

Returns `kind` – kind of recording in [monitorTrend, monitorWave, taphTrend, telvet]

Return type *str*

`anesplot.record_main.select_wave(waves, num=1)`

select the recording type:

Returns `kind` – kind of recording in [monitorTrend, monitorWave, taphTrend, telvet]

Return type *str*

`anesplot.record_main.build_param_dico(file=None, asource=None, pathdico={'cwd':`

`'/Users/cdesbois/pg/chrisPg/anesthPlot', 'data':`

`'/Users/cdesbois/enva/clinique/recordings/anesthRecords/onPanelPcRecorded',`

`'recordMain': '/Users/cdesbois/pg/chrisPg/anesthPlot/anesplot',`

`'root': '/Users/cdesbois', 'sBg': '/Users/cdesbois/ownCloud', 'sFig':`

`'/Users/cdesbois/ownCloud', 'save': '/Users/cdesbois/ownCloud',`

`'utils': '/Users/cdesbois/pg'})`

initialise a dict save parameters —> TODO see min vs sec

Parameters

- **file** (*str*) – the recording filename

- **source** (*str*) – the origin of the recording

Returns

dico –

a dictionary describing the situation [item, xmin, xmax, ymin, ymax, path, unit, save, memo, file, source]

Return type dict

`anesplot.record_main.plot_trenddata(df, header, param_dico)`

clinical main plots of a trend recordings

parameters df : pdDataframe

recorded data (MonitorTrend.data)

header [dict] recording parameters (MonitorTrend.header)

param_dico [dict] plotting parameters (MonitorTrend.param)

Returns **afig_dico**

Return type dict of name:fig

class `anesplot.record_main.Waves(filename=None)`

Bases: object

the base object to store the records.

class `anesplot.record_main.SlowWave(filename=None)`

Bases: [anesplot.record_main.Waves](#)

class for slowWaves = trends

file [str] short name

filename [str] long name

clean_trend : **external**

clean the data

show_graphs : **external**

plot clinical main plots

clean_trend()

clean the data, remove irrelevant, input = self.data, output = pandas dataframe nb doesnt change the obj.data in place

show_graphs()

basic clinical plots

class `anesplot.record_main.MonitorTrend(filename=None, load=True)`

Bases: [anesplot.record_main.SlowWave](#)

monitor trends recordings:

input = filename : path to file load = boolean to load data (default is True)

file [str] short name

filename [str] long name

header [dict] record parameters

source [str] recording apparatus (default = 'monitor')

fs [float] sampling rate

param [dict] display parameters

clean_trend [external] clean the data

show_graphs [external] plot clinical main plots

class `anesplot.record_main.TaphTrend(filename=None)`

Bases: `anesplot.record_main.SlowWave`

taphonius trends recordings

input ... FILLME

attributes ... FILLME

load_header()

load the header -> pandas.dataframe

extract_taph_events()

extract Taph events

Parameters `data` (`pandas dataframe`) – record df form taphonius recording)

Returns events dataframe

Return type eventdf pandas dataframe

class `anesplot.record_main.FastWave(filename=None)`

Bases: `anesplot.record_main.Waves`

class for Fastwaves = continuous recordings.

plot_wave(tracesList=None)

simple choose and plot for a wave input = none -> GUI, or list of waves to plot (max=2)

define_a_roi()

define a ROI.

class `anesplot.record_main.TeleVetWave(filename=None)`

Bases: `anesplot.record_main.FastWave`

class to organise teleVet recordings transformed to csv files.

class `anesplot.record_main.MonitorWave(filename=None, load=True)`

Bases: `anesplot.record_main.FastWave`

class to organise monitorWave recordings. input : filename = path to file load = boolean to load data (default is True)

attibutes ... FILLME

methods ... FILLME

`anesplot.record_main.main()`

MODULES

3.1 anesplot package

3.1.1 Subpackages

anesplot.config package

Submodules

anesplot.config.build_recordRc module

build a 'recordRc.yaml' configuration file to adapt to a specific computer location at the root of anesplot

- input <-> 'data' : to load the records
- output <-> 'save' : to save the plots

```
anesplot.config.build_recordRc.filedialog(kind="", direc-  
                                         tory='/Users/cdesbois/pg/chrisPg/anesthPlot/anesplot/config',  
                                         for_open=True, fmt="", is_folder=False)  
    general dialog function.
```

```
anesplot.config.build_recordRc.read_config()  
    locate & load the yaml file.
```

```
anesplot.config.build_recordRc.write_configfile(path)  
    record the yaml file.
```

```
anesplot.config.build_recordRc.main()  
    main function for script execution.
```

anesplot.config.load_recordRc module

load an already generated 'recordRc.yaml' configuration file

- input <-> 'data' : to load the records
- output <-> 'save' : to save the plots

`anesplot.config.load_recordRc.build_paths()`

read the yaml configuration file.

`anesplot.config.load_recordRc.adapt_with_syspath(path_dico)`

add the folder location to the system path.

Module contents

anesplot.loadrec package

Submodules

anesplot.loadrec.explore module

Created on Thu Mar 12 16:52:13 2020

@author: cdesbois

`anesplot.loadrec.explore.gui_choosefile(paths=None)`

select a file via a dialog and return the file name.

anesplot.loadrec.loadmonitor_trendrecord module

Created on Wed Jul 24 13:43:26 2019 @author: cdesbois

load a monitor trend recording:

- choose a file
- load the header to a dictionary
- load the data into a pandas dataframe

`anesplot.loadrec.loadmonitor_trendrecord.choosefile_gui(dir_path=None)`

select a file using a dialog.

Parameters `dir_path` (*str*) – optional location of the data (`paths['data']`)

Returns filename (full path)

Return type str

`anesplot.loadrec.loadmonitor_trendrecord.loadmonitor_trendheader(filename)`

load the file header.

Parameters `filename` (*str*) – full name of the file

Returns header

Return type dict

`anesplot.loadrec.loadmonitor_trendrecord.loadmonitor_trenddata(filename, header)`

load the monitor trend data

Parameters

- `filename` (*str*) – filename

- **header** (*dict*) – fileheader

Returns df = trends data

Return type pandas.DataFrame

anesplot.loadrec.loadmonitor_waverecord module

Created on Wed Jul 24 14:56:58 2019 @author: cdesbois

load a monitor wave recording:

- choose a file
- load the header to a pandas dataframe
- load the date into a pandas dataframe

anesplot.loadrec.loadmonitor_waverecord.**choosefile_gui**(*dir_path=None*)
select a file using a dialog.

Parameters **dir_path** (*str*) – optional location of the data (paths['data'])

Returns filename (full path)

Return type str

anesplot.loadrec.loadmonitor_waverecord.**loadmonitor_waveheader**(*filename*)
load the wave file header.

Parameters **filename** (*str*) – full name of the file

Returns header

Return type pandas.DataFrame

anesplot.loadrec.loadmonitor_waverecord.**loadmonitor_wavedata**(*filename*)
load the monitor wave csvDataFile.

Parameters **filename** (*str*) – full name of the file

Returns df = trends data

Return type pandas.DataFrame

anesplot.loadrec.loadtaph_trendrecord module

Created on Wed Jul 24 15:30:07 2019 @author: cdesbois

load a taphonius data recording:

- choose a file
- load the patient datafile to a dictionary
- load the physiological date into a pandas dataframe

anesplot.loadrec.loadtaph_trendrecord.**choosefile_gui**(*dir_path=None*)
select a file using a dialog.

Parameters `dir_path` (*str*) – optional location of the data (paths[‘data’])

Returns filename (full path)

Return type str

`anesplot.loadrec.loadtaph_trendrecord.loadtaph_trenddata(filename)`

load the taphoniusData trends data.

Parameters `filename` (*str*) – fullname

Returns df = trends data

Return type pandas.DataFrame

`anesplot.loadrec.loadtaph_trendrecord.loadtaph_patientfile(headername)`

load the taphonius patient.csv file

Parameters `headername` (*str*) – fullname

Returns descr = patient_data

Return type dict

anesplot.loadrec.loadtelevet module

Created on Wed Jul 31 16:22:06 2019 @author: cdesbois

load televet exported (csv) data: to be developped

`anesplot.loadrec.loadtelevet.choosefile_gui(dir_path=None)`

select a file using a dialog.

Parameters `dir_path` (*str*) – optional location of the data (paths[‘data’])

Returns filename (full path)

Return type str

`anesplot.loadrec.loadtelevet.loadtelevet(file=None, all_traces=False)`

load the televetCsvExportedFile.

Parameters

- **file** (*str*) – name of the file
- **all_traces** (*bool*) – load all the derivations

Returns df = recorded traces

Return type pandas.DataFrame

Module contents

anesplot.plot package

Submodules

anesplot.plot.trend_plot module

Created on Tue Apr 19 09:08:56 2016 @author: cdesbois

collection of functions to plot the trend data

`anesplot.plot.trend_plot.color_axis(ax, spine='bottom', color='r')`
change the color of the label & tick & spine.

Parameters

- **ax** (`matplotlib.pyplot.axis`) – the axis
- **spine** (`str`) – optional location in ['bottom', 'left', 'top', 'right']
- **colors** (`str`) – optional color

`anesplot.plot.trend_plot.append_loc_to_fig(ax, dt_list, label='g')`
append vertical lines to indicate a location 'eg: arterial blood gas'

Parameters

- **ax** (`matplotlib.pyplot.axis`) – the axis
- **dt_list** (`[datetime]`) – list of datetime values
- **label** (`str`) – a key to add to the label (default is 'g')

Returns `res` a dictionary containing the locations

Return type `dict`

`anesplot.plot.trend_plot.plot_header(descr, param={'save': False})`
plot the header of the file.

Parameters

- **descr** (`dict`) – header of the recording
- **param** (`dict`) – dictionary of parameters

Returns `fig` plot of the header

Return type `pyplot.figure`

`anesplot.plot.trend_plot.hist_cardio(data, param={})`
mean arterial pressure histogramme using matplotlib.

Parameters

- **data** (`pandas.DataFrame`) – the recorded trends data (keys used : 'ip1m' and 'hr'),
- **param** (`dict`) – parameters (save=boolean, 'path': path to directory)

Returns `fig` matplotlib.pyplot.figure

`anesplot.plot.trend_plot.plot_one_over_time(x, y, colour)`
plot y over x using colour

`anesplot.plot.trend_plot.hist_co2_iso(data, param={})`

CO2 and iso histogramme (NB CO2 should have been converted from % to mmHg)

Parameters

- **data** (*pandas.DataFrame*) – the trends recorded data
- **param** (*dict*) – dictionary of parameters

Returns fig pyplot.figure

`anesplot.plot.trend_plot.cardiovasc(data, param={})`

cardiovascular plot

Parameters

- **data** (*pandas.DataFrame*) – the recorded trends data keys used :['ip1s', 'ip1m', 'ip1d', 'hr']
- **param** (*dict*) – dict(save: boolean, path['save'], xmin, xmax, unit, dtime = boolean for time display in HH:MM format)

Returns fig= pyplot.figure

`anesplot.plot.trend_plot.cardiovasc_p1p2(data, param={})`

cardiovascular plot with central venous pressure (p2)

Parameters

- **data** (*pandas.DataFrame*) – the trends recorded data keys used :['ip1s', 'ip1m', 'ip1d', 'hr', 'ip2s', 'ip2m', 'ip2d']
- **param** (*dict*) – dict(save: boolean, path['save'], xmin, xmax, unit, dtime = boolean for time display in HH:MM format)

Returns fig= pyplot.figure

`anesplot.plot.trend_plot.co2iso(data, param={})`

anesth plot (CO2/iso)

Parameters

- **data** (*pandas.DataFrame*) – the recorded data keys used :['ip1s', 'ip1m', 'ip1d', 'hr']
- **param** (*dictionary*) – dict(save: boolean, path['save'], xmin, xmax, unit, dtime = boolean for time display in HH:MM format)

:returns fig= pyplot.figure

`anesplot.plot.trend_plot.func(ax, x, y1, y2, color='tab:blue', x0=38)`

`anesplot.plot.trend_plot.co2o2(data, param)`

respiratory plot (CO2 and Iso)

Parameters

- **data** (*pandas.DataFrame*) – recorded trends data keys used :['ip1s', 'ip1m', 'ip1d', 'hr']
- **param** (*dict*) – dict(save: boolean, path['save'], xmin, xmax, unit, dtime = boolean for time display in HH:MM format)

Returns fig= pyplot.figure

`anesplot.plot.trend_plot.ventil(data, param)`

plot ventilation parameters (.tvInsp, .pPeak, .pPlat, .peep, .minVexp, .co2RR, .co2exp)

Parameters

- **data** (*pandas.DataFrame*) – recorded data, keys used :['ip1s', 'ip1m', 'ip1d', 'hr']
- **param** (*dict*) – dict(save: boolean, path['save'], xmin, xmax, unit, dtime = boolean for time display in HH:MM format)

Returns fig= pyplot.figure

anesplot.plot.trend_plot.**recrut**(*data, param*)
display a recrut manoeuver (.pPeak, .pPlat, .peep, .tvInsp)

Parameters

- **data** (*pandas.DataFrame*) – recorded data keys used :['ip1s', 'ip1m', 'ip1d', 'hr']
- **param** (*dict*) – dict(save: boolean, path['save'], xmin, xmax, unit, dtime = boolean for time display in HH:MM format)

:returns fig= pyplot.figure

anesplot.plot.trend_plot.**ventil_cardio**(*data, param*)
build ventilation and cardiovascular plot

Parameters

- **data** (*pandas.DataFrame*) – teh recorded trends data keys used :['ip1s', 'ip1m', 'ip1d', 'hr']
- **param** (*dict*) – dict(save: boolean, path['save'], xmin, xmax, unit, dtime = boolean for time display in HH:MM format)

Returns fig= pyplot.figure

anesplot.plot.trend_plot.**save_distri**(*data, path*)
save as 'O_.' the 4 distributions graphs for cardiovasc annd respi

anesplot.plot.trend_plot.**fig_memo**(*path, fig_name*)
append latex citation commands in a txt file inside the fig folder create the file iif it doesn't exist

anesplot.plot.wave_plot module

Created on Tue Apr 19 09:08:56 2016

@author: cdesbois

anesplot.plot.wave_plot.**color_axis**(*ax, spine='bottom', color='r'*)
change the color of the label & tick & spine.

Parameters

- **ax** (*matplotlib.pyplot.axis*) – the axis
- **spine** (*str*) – optional location in ['bottom', 'left', 'top', 'right']
- **colors** (*str*) – optional color

anesplot.plot.wave_plot.**plot_wave**(*data, keys=[], param={}*)
plot the waves recorded (from as5)

Parameters

- **data** (*pandas.DataFrame*) – the recorded trends data
- **keys** (*list*) – one or two in ['wekg', 'ECG', 'wco2', 'wawp', 'wflow', 'wap']
- **{mini}** (*dict*) – limits in point value (index), maxi: limits in point value (index)}

Returns fig plt.figure the plot

Returns lines plt.line2D the line to animate

(Nb plot data/index, but the xscale is indicated as sec)

Module contents

Created on Tue Apr 19 09:08:56 2016

functions to plot the trend data

@author: cdesbois

anesplot.treatrec package

Submodules

anesplot.treatrec.build_anim module

anesplot.treatrec.clean_data module

Created on Wed Jul 31 16:05:29 2019

@author: cdesbois

`anesplot.treatrec.clean_data.clean_trenddata(df)`
remove artifacts in the recorded trends

anesplot.treatrec.ekg_to_hr module

Created on Wed Feb 12 16:52:00 2020 @author: cdesbois

function used to treat an EKG signal and extract the heart rate typically (copy, paste and execute line by line)

0. after

```
import anesplot.record_main as rec
from treatrec import ekg_to_hr as tohr
```

1. load the data in a pandas dataframe:

(through classes rec.MonitorTrend & rec.MonitorWave)

```
trendname = '' # fullname
or
trendname = rec.choosefile_gui()
```

```
wavename = rec.trendname_to_wavename(trendname)
-
# load the data
trends = rec.MonitorTrend(trendname)
waves = rec.MonitorWave(wavename)
-
# format the name
name = trends.header['Patient Name'].title().replace(' ', '')
name = name[0].lower() + name[1:]
```

2. treat the ekg wave:

- get parameters
- build a dataframe to work with (waves)
- low pass filtering
- build the beat locations (beat based dataframe):

```
params = waves.param
ekg_df = pd.DataFrame(waves.data.wekg)
ekg_df['wekg_lowpass'] = rec.wf.fix_baseline_wander(ekg_df.wekg,
                                                    waves.param['fs'])
beat_df = tohr.detect_beats(ekg_df.wekg_lowpass, mult=1)
```

3. perform the manual adjustments required:

- based on a graphical display of beat locations, an rr values
- build a container for the manual corrections:

```
figure = tohr.plot_beats(ekg_df.wekg_lowpass, beat_df)
to_change_df = pd.DataFrame(columns=beat_df.columns.insert(0, 'action'))
```

- remove or add peaks : zoom on the figure to observe only one peak, then:

```
to_change_df = tohr.remove_beat(beat_df, ekg_df, to_change_df, figure)
or
to_change_df = tohr.append_beat(beat_df, ekg_df, to_change_df, figure,
                                yscale=1)
```

- combine to update the beat_df with the manual changes:

```
beat_df = tohr.update_beat_df(beat_df, to_change_df,
                              path_to_file="", from_file=False)
```

- save the peaks locations:

```
tohr.save_beats(beat_df, to_change_df, savename='', savepath=None)
(# or reload
beat_df = pd.read_hdf('beatDf.hdf', key='beatDf') )
```

4. go from points values to continuous time:

```
beat_df = tohr.compute_rr(beat_df)
ahr_df = tohr.interpolate_rr(beat_df)
tohr.plot_rr(ahr_df, params)
```

5. append intantaneous heart rate to the initial data:

```
ekg_df = tohr.append_rr_and_ihr_to_wave(ekg_df, ahr_df)
waves.data = tohr.append_rr_and_ihr_to_wave(waves.data, ahr_df)
trends.data = tohr.append_ihr_to_trend(trends.data, waves.data, ekg_df)
```

6. save:

```
tohr.save_trends_data(trends.data, savename=name, savepath='data')
tohr.save_waves_data(waves.data, savename=name, savepath='data')
```

```
anesplot.treatrec.ekg_to_hr.detect_beats(ser, fs=300, species='horse', mult=1)
    detect the peak locations
```

Parameters

- **ser** (*pandas.series*) – the data
- **fs** (*integer*) – sampling frequency
- **species** (*string*) – in [horse]
- **mult** (*float*) – correction / 1 for qRs amplitude

Returns df=pandas.DataFrame

```
anesplot.treatrec.ekg_to_hr.plot_beats(ecg, beats)
    plot ecg waveform + beat location
```

```
anesplot.treatrec.ekg_to_hr.append_beat(beatdf, ekgdf, tochange_df, fig, lim=None, yscale=1)
    locate the beat in the figure, append to a dataframe['toAppend']
```

Parameters

- **beatdf** (*pandas.DataFrame*) – contains the point based location (pLocs)
- **dataframe ekgdf** (*pandas*) – contains the wave recording ((wekg_lowpass)
- **tochange_df** (*pandas.DataFrame*) – to store the beats toAppend or toRemove
- **fig** (*pyplot.Figure*) – figure to find time limits
- **lim** (*integer*) – ptBasedLim optional to give it manually
- **yscale** (*float*) – amplitude mutliplication factor for detection (default=1)

Returns tochange_df: incremented changedf (pt location)

Return type pandasDataFrame

methods :

locate the beat in the figure, append to a dataframe['toAppend'] 0.: if not present : build a dataframe:

```
>>> to_change_df = pd.DataFrame(columns=['toAppend', 'toRemove'])
```

1.: locate the extra beat in the figure (cf plot_beats()) and zoom to observe only a negative peak

2.: call the function:

```
>>> to_change_df = remove_beat(beatdf, ekgdf, tochange_df, fig)
-> the beat parameters will be added the dataframe
```

.in the end of the manual check, update the beat_df

- first : save beat_df and to_change_df
- **second** [run:]

```
>>> beat_df = update_beat_df()
```

`anesplot.treatrec.ekg_to_hr.remove_beat(beatdf, ekgdf, tochange_df, fig, lim=None)`
locate the beat in the figure, append to a dataframe['toRemove']

0.: if not present build a dataframe:

```
>>> to_change_df = pd.DataFrame(columns=['toAppend', 'toRemove'])
```

1.: locate the extra beat in the figure (cf plot_beats()) and zoom to observe only a negative peak

2.: call the function:::

```
>>> to_change_df = remove_beat(beatdf, ekgdf, tochange_df, fig)
-> the beat parameters will be added the dataframe
```

.(in the end of the manual check, update the beat_df

- first : save beat_df and to_change_df
- **second** [run]

```
>>> beat_df = update_beat_df()
```

`anesplot.treatrec.ekg_to_hr.save_beats(beatdf, tochangedf, savename="", savepath=None)`
save the beats locations as csv and hd5 file

Parameters

- **beatde** (*pd.dataframes*) –
- **savepath** (*path to save in*) –

`anesplot.treatrec.ekg_to_hr.update_beat_df(beatdf, tochangedf, path_to_file="", from_file=False)`
implement in the beat location the manual corrections fromFile = True force the disk loading of the dataframes

`anesplot.treatrec.ekg_to_hr.compute_rr(beatdf, fs=None)`
compute rr intervals (from pt to time)

Parameters

- **beatdf** (*pd.DataFrame*) – with 'pLoc'
- **fs** (*integer*) – sampling frequency

Returns with: 'rr' = rr duration 'rrDiff' = rrVariation 'rrSqDiff' = rrVariation^2

Return type pd.DataFrame

`anesplot.treatrec.ekg_to_hr.interpolate_rr(beatdf, kind=None)`
 interpolate the beat_df (pt -> time values)

Parameters

- **beatDf** (*pd.DataFrame*) –
- **kind** (*str*) – 'linear' or 'cubic'(default)

Returns 'espts' = evenly spaced points 'rrInterpol' = interpolated rr

Return type pdDataframe with evenly spaced data

`anesplot.treatrec.ekg_to_hr.plot_rr(ahr_df, param, HR=False)`
 plot RR vs pt values + rrSqDiff

Parameters

- **pdDataFrame** (*hr_df*) –
- **params** – dict containing 'fs' as key

`anesplot.treatrec.ekg_to_hr.append_rr_and_ihr_to_wave(wave, ahrdf)`
 append rr and ihr to the waves based on pt value (ie index)

`anesplot.treatrec.ekg_to_hr.plot_agreement(trenddf)`
 plot ip1HR & ihr to check agreement

`anesplot.treatrec.ekg_to_hr.append_ihr_to_trend(trenddf, wavedf, ekgdf)`
 append 'ihr' (instantaneous heart rate) to the trends

`anesplot.treatrec.ekg_to_hr.save_trends_data(trenddf, savename="", savepath=None)`
 save the trends data to a csv and hd5 file, including an ihr column

Parameters

- **trenddf** (*pd.dataframes*) –
- **savename** (*str*) –
- **savepath** (*str*) – path to save in (default= current working directory)

`anesplot.treatrec.ekg_to_hr.save_waves_data(wavedf, savename="", savepath=None)`
 save the trends data to a csv and hd5 file, including an ihr column

Parameters

- **trenddf** (*pd.dataframes*) –
- **savename** (*str*) – savepath : path to save in (default= current working directory)

anesplot.treatrec.extract_hypotension module

Spyder Editor

This is a temporary script file.

`anesplot.treatrec.extract_hypotension.extract_hypotension(trends, pamin=70)`
return a dataframe with the beginning and ending phses of hypotension

Parameters

- **trends** (*MonitorTrend object*) –
- **pamin** (*float= threshold de define hypotension on mean arterial pressure*) –
- **is 70** (*((default)*) –

Returns **durdf** – transitions (up and down, in seconds from beginning) and duration in the hypotension state (in seconds)

Return type pandas DataFrame containing

`anesplot.treatrec.extract_hypotension.plot_hypotension(trends, durdf, durmin=15, pamin=70)`
plot the hupotentions phases

Parameters

- **trends** (*TYPE*) – DESCRIPTION.
- **durdf** (*TYPE*) – DESCRIPTION.
- **durmin** (*TYPE, optional*) – DESCRIPTION. The default is 15.

Returns **fig** – DESCRIPTION.

Return type TYPE

`anesplot.treatrec.extract_hypotension.scatter_length_meanhypo(trends, durdf)`
draw a scatter plot (hypotensive arterial value vs duration of hypotension) :param trends: :type trends: MonitorTrend :param durdf: :type durdf: pandas dataframe containing the value and duration

Returns **fig**

Return type matplotlib.pyplot figure

`anesplot.treatrec.extract_hypotension.plot_all_dir_hypo(dirname=None, scatter=False)`
walk throught the folder and plot the values

anesplot.treatrec.hr_to_hrv module

`anesplot.treatrec.hr_to_hrv.build_hrv_limits(spec='horse')`
return a dico containing HRV limits (VLF, LF, HF) input : spec in ['horse', 'man']

anesplot.treatrec.wave_func module

Created on Fri Dec 8 12:46:41 2017

@author: cdesbois

`anesplot.treatrec.wave_func.fix_baseline_wander(data, fs=500)`

BaselineWanderRemovalMedian.m from ecg-kit. Given a list of amplitude values (data) and sample rate (sr), it applies two median filters to data to compute the baseline. The returned result is the original data minus this computed baseline.

`anesplot.treatrec.wave_func.rol_mean(ser, win_lengh=1, fs=500)`

returns a rolling mean of a RR serie

Parameters

- **pd.Serie** (*ser*) –
- **win_lengh** (*integer*) – window lenght for averaging (in sec),
- **fs** (*int*) – sampling frequency

`anesplot.treatrec.wave_func.return_points(df, fig)`

return a tuple containing the point values of ROI

Parameters

- **df** (*anesthesia record dataframe*) –
- **fig** (*pyplot.figure*) –

Returns ROI

Return type dict

`anesplot.treatrec.wave_func.restrict_time_area(df1, mini=None, maxi=None)`

return a new dataframe with reindexation

Parameters

- **df1** (*pandas.DataFrame*) –
- **mini** (*integer*) – miniPointValue
- **maxi** (*integer*) – maxiPointValue

Returns

Return type pandas.DataFrame

Module contents

3.1.2 Submodules

3.1.3 Module contents

INDICES AND TABLES

- `genindex`
- `modindex`
- `search`

PYTHON MODULE INDEX

a

- `anesplot`, 20
- `anesplot.config`, 8
 - `anesplot.config.build_recordRc`, 7
 - `anesplot.config.load_recordRc`, 7
- `anesplot.loadrec`, 11
 - `anesplot.loadrec.explore`, 8
 - `anesplot.loadrec.loadmonitor_trendrecord`, 8
 - `anesplot.loadrec.loadmonitor_waverecord`, 9
 - `anesplot.loadrec.loadtaph_trendrecord`, 9
 - `anesplot.loadrec.loadtelevet`, 10
- `anesplot.plot`, 14
 - `anesplot.plot.trend_plot`, 11
 - `anesplot.plot.wave_plot`, 13
- `anesplot.record_main`, 3
- `anesplot.treatrec`, 20
 - `anesplot.treatrec.clean_data`, 14
 - `anesplot.treatrec.ekg_to_hr`, 14
 - `anesplot.treatrec.extract_hypotension`, 19
 - `anesplot.treatrec.hr_to_hrv`, 19
 - `anesplot.treatrec.wave_func`, 20

A

`adapt_with_syspath()` (in module *anesplot.config.load_recordRc*), 8
`anesplot`
 module, 20
`anesplot.config`
 module, 8
`anesplot.config.build_recordRc`
 module, 7
`anesplot.config.load_recordRc`
 module, 7
`anesplot.loadrec`
 module, 11
`anesplot.loadrec.explore`
 module, 8
`anesplot.loadrec.loadmonitor_trendrecord`
 module, 8
`anesplot.loadrec.loadmonitor_waverecord`
 module, 9
`anesplot.loadrec.loadtaph_trendrecord`
 module, 9
`anesplot.loadrec.loadtelevet`
 module, 10
`anesplot.plot`
 module, 14
`anesplot.plot.trend_plot`
 module, 11
`anesplot.plot.wave_plot`
 module, 13
`anesplot.record_main`
 module, 3
`anesplot.treatrec`
 module, 20
`anesplot.treatrec.clean_data`
 module, 14
`anesplot.treatrec.ekg_to_hr`
 module, 14
`anesplot.treatrec.extract_hypotension`
 module, 19
`anesplot.treatrec.hr_to_hrv`
 module, 19
`anesplot.treatrec.wave_func`

 module, 20
`append_beat()` (in module *anesplot.treatrec.ekg_to_hr*), 16
`append_ihr_to_trend()` (in module *anesplot.treatrec.ekg_to_hr*), 18
`append_loc_to_fig()` (in module *anesplot.plot.trend_plot*), 11
`append_rr_and_ihr_to_wave()` (in module *anesplot.treatrec.ekg_to_hr*), 18

B

`build_hrv_limits()` (in module *anesplot.treatrec.hr_to_hrv*), 19
`build_param_dico()` (in module *anesplot.record_main*), 3
`build_paths()` (in module *anesplot.config.load_recordRc*), 7

C

`cardiovasc()` (in module *anesplot.plot.trend_plot*), 12
`cardiovasc_p1p2()` (in module *anesplot.plot.trend_plot*), 12
`check()` (in module *anesplot.record_main*), 4
`choosefile_gui()` (in module *anesplot.loadrec.loadmonitor_trendrecord*), 8
`choosefile_gui()` (in module *anesplot.loadrec.loadmonitor_waverecord*), 9
`choosefile_gui()` (in module *anesplot.loadrec.loadtaph_trendrecord*), 9
`choosefile_gui()` (in module *anesplot.loadrec.loadtelevet*), 10
`choosefile_gui()` (in module *anesplot.record_main*), 3
`clean_trend()` (*anesplot.record_main.SlowWave* method), 4
`clean_trenddata()` (in module *anesplot.treatrec.clean_data*), 14
`co2iso()` (in module *anesplot.plot.trend_plot*), 12
`co2o2()` (in module *anesplot.plot.trend_plot*), 12
`color_axis()` (in module *anesplot.plot.trend_plot*), 11
`color_axis()` (in module *anesplot.plot.wave_plot*), 13
`compute_rr()` (in module *anesplot.treatrec.ekg_to_hr*), 17

D

`define_a_roi()` (*anesplot.record_main.FastWave method*), 5
`detect_beats()` (*in module anesplot.treatrec.ekg_to_hr*), 16

E

`extract_hypotension()` (*in module anesplot.treatrec.extract_hypotension*), 19
`extract_taph_events()` (*anesplot.record_main.TaphTrend method*), 5

F

`FastWave` (*class in anesplot.record_main*), 5
`fig_memo()` (*in module anesplot.plot.trend_plot*), 13
`filedialog()` (*in module anesplot.config.build_recordRc*), 7
`fix_baseline_wander()` (*in module anesplot.treatrec.wave_func*), 20
`func()` (*in module anesplot.plot.trend_plot*), 12

G

`gui_choosefile()` (*in module anesplot.loadrec.explore*), 8

H

`hist_cardio()` (*in module anesplot.plot.trend_plot*), 11
`hist_co2_iso()` (*in module anesplot.plot.trend_plot*), 12

I

`interpolate_rr()` (*in module anesplot.treatrec.ekg_to_hr*), 18

L

`load_header()` (*anesplot.record_main.TaphTrend method*), 5
`loadmonitor_trenddata()` (*in module anesplot.loadrec.loadmonitor_trendrecord*), 8
`loadmonitor_trendheader()` (*in module anesplot.loadrec.loadmonitor_trendrecord*), 8
`loadmonitor_wavedata()` (*in module anesplot.loadrec.loadmonitor_waverecord*), 9
`loadmonitor_waveheader()` (*in module anesplot.loadrec.loadmonitor_waverecord*), 9
`loadtaph_patientfile()` (*in module anesplot.loadrec.loadtaph_trendrecord*), 10
`loadtaph_trenddata()` (*in module anesplot.loadrec.loadtaph_trendrecord*), 10
`loadtelevet()` (*in module anesplot.loadrec.loadtelevet*), 10

M

`main()` (*in module anesplot.config.build_recordRc*), 7
`main()` (*in module anesplot.record_main*), 5
module
 `anesplot`, 20
 `anesplot.config`, 8
 `anesplot.config.build_recordRc`, 7
 `anesplot.config.load_recordRc`, 7
 `anesplot.loadrec`, 11
 `anesplot.loadrec.explore`, 8
 `anesplot.loadrec.loadmonitor_trendrecord`, 8
 `anesplot.loadrec.loadmonitor_waverecord`, 9
 `anesplot.loadrec.loadtaph_trendrecord`, 9
 `anesplot.loadrec.loadtelevet`, 10
 `anesplot.plot`, 14
 `anesplot.plot.trend_plot`, 11
 `anesplot.plot.wave_plot`, 13
 `anesplot.record_main`, 3
 `anesplot.treatrec`, 20
 `anesplot.treatrec.clean_data`, 14
 `anesplot.treatrec.ekg_to_hr`, 14
 `anesplot.treatrec.extract_hypotension`, 19
 `anesplot.treatrec.hr_to_hrv`, 19
 `anesplot.treatrec.wave_func`, 20
`MonitorTrend` (*class in anesplot.record_main*), 4
`MonitorWave` (*class in anesplot.record_main*), 5

P

`plot_agreement()` (*in module anesplot.treatrec.ekg_to_hr*), 18
`plot_all_dir_hypo()` (*in module anesplot.treatrec.extract_hypotension*), 19
`plot_beats()` (*in module anesplot.treatrec.ekg_to_hr*), 16
`plot_header()` (*in module anesplot.plot.trend_plot*), 11
`plot_hypotension()` (*in module anesplot.treatrec.extract_hypotension*), 19
`plot_monitorwave_data()` (*in module anesplot.record_main*), 4
`plot_one_over_time()` (*in module anesplot.plot.trend_plot*), 11
`plot_rr()` (*in module anesplot.treatrec.ekg_to_hr*), 18
`plot_trenddata()` (*in module anesplot.record_main*), 4
`plot_wave()` (*anesplot.record_main.FastWave method*), 5
`plot_wave()` (*in module anesplot.plot.wave_plot*), 13

R

`read_config()` (*in module anesplot.config.build_recordRc*), 7
`recrut()` (*in module anesplot.plot.trend_plot*), 13

remove_beat() (in module *anesplot.treatrec.ekg_to_hr*), 17
 restrict_time_area() (in module *anesplot.treatrec.wave_func*), 20
 return_points() (in module *anesplot.treatrec.wave_func*), 20
 rol_mean() (in module *anesplot.treatrec.wave_func*), 20

S

save_beats() (in module *anesplot.treatrec.ekg_to_hr*), 17
 save_distri() (in module *anesplot.plot.trend_plot*), 13
 save_trends_data() (in module *anesplot.treatrec.ekg_to_hr*), 18
 save_waves_data() (in module *anesplot.treatrec.ekg_to_hr*), 18
 scatter_length_meanhypo() (in module *anesplot.treatrec.extract_hypotension*), 19
 select_type() (in module *anesplot.record_main*), 3
 select_wave() (in module *anesplot.record_main*), 3
 show_graphs() (*anesplot.record_main.SlowWave* method), 4
 SlowWave (class in *anesplot.record_main*), 4

T

TaphTrend (class in *anesplot.record_main*), 5
 TelevetWave (class in *anesplot.record_main*), 5
 trendname_to_wavename() (in module *anesplot.record_main*), 3

U

update_beat_df() (in module *anesplot.treatrec.ekg_to_hr*), 17

V

ventil() (in module *anesplot.plot.trend_plot*), 12
 ventil_cardio() (in module *anesplot.plot.trend_plot*), 13

W

Waves (class in *anesplot.record_main*), 4
 write_configfile() (in module *anesplot.config.build_recordRc*), 7