

---

# **anesthPlot**

***Release beta***

**Christophe Desbois**

**Jan 26, 2022**



## **MAIN<sub>s</sub>CRIP*T* :**

<b>1</b>	<b>Welcome to anesthPlots documentation!</b>	<b>1</b>
<b>2</b>	<b>main script</b>	<b>3</b>
<b>3</b>	<b>modules</b>	<b>5</b>
<b>4</b>	<b>Indices and tables</b>	<b>21</b>
	<b>Python Module Index</b>	<b>23</b>
	<b>Index</b>	<b>25</b>



## WELCOME TO ANESTHPLOTS DOCUMENTATION!

anesthPlot is a python package developped to extract, manipulate and plots anesthesia data recorded from the Monitor Software to be used mostly in a teaching environment.

**Warning:** This project is:

- a work in progres
- the processes are mainly focused on horses anesthesia
- in our environment the data recorded came from either
  - an as3 or as5 anesthesia monitor (ekg, invasive pressure, etCO2, halogenate, spirometry)
  - a Taphonius equine ventilator
  - (some ekg data extracted using a Televet holter system)

### 1.1 Features

- you can **load** recordings from a trend or a wave file
  - **from command line:**

```
python anesthPlot/anesplot/__main__.py
-> will open an GUI choose menu to select the recording
(monitorTrend, taphoniusTrend, monitorWave, televetWave(export))
```

- \* will build a **standard debriefing** (trends) **plot series** (script usage)
  - global histograms (cardiovascular and anesthesia summary)
  - cardiovascular trends time based plots
  - respiratory trends time based plots
  - anesthesia trends time based plots
- \* or will build a **plot for wave** recording
  - one or two waves on the same plot (script usage)
- you can also use this code as a **python package**
  - **usage :**

```
import anesplot.record_main as rec
trendname = 'a_full_path_to_csv_file'
# nb if no filename is provided, a chooseFile Gui will be called to choose_
↪ the file
trends = rec.MonitorTrend(trendname)
#(you can also use trends = rec.taphTrend())
wavename = rec.trendname_to_wavename(trendname)
waves = rec.MonitorWave(trends)

trends.show_graphs() # -> set of plots for 'clinical' debriefing purposes

waves.plot_waves() # -> one or two traces
# ... adjust manually the scales of the display
waves.define_a_roi() # -> to register the plotting scales
waves.animate_fig() #-> to build an animation using these parameters
```

- additional functions are available to extract instaneous heart rate
  - \* see anesplot/treatrec/ekg\_to\_hr.py

## MAIN SCRIPT

### 2.1 anesplot.record\_main module

main script/module to load and display an anesthesia record

**can be runned as a script::** python record\_main.py

**or imported as a package::** import anesplot.record\_main as rec %gui qt5 (required only to use the dialogs if using spyder) trends = rec.MonitorTrend() waves = rec.MonitorWave(rec.trendname\_to\_wavename(trends.filename))

anesplot.record\_main.**choosefile\_gui**(dirname: *Optional[str] = None*) → str  
Select a file via a dialog and return the (full) filename.

**Parameters** **dir\_path** (str) – location to place the gui (generally paths[data]) else home

**Returns** **fname[0]** – filename

**Return type** str

anesplot.record\_main.**trendname\_to\_wavename**(name: str) → str  
just compute the supposed name

anesplot.record\_main.**select\_type**(question: *Optional[str] = None*, items: *Optional[list] = None*, num: *int = 0*) → str

select the recording type:

**Returns** **kind** – kind of recording in [monitorTrend, monitorWave, taphTrend, telvet]

**Return type** str

anesplot.record\_main.**select\_wave\_to\_plot**(waves: list, num=1) → str  
select the recording type:

**Returns** **kind** – kind of recording in [monitorTrend, monitorWave, taphTrend, telvet]

**Return type** str

anesplot.record\_main.**plot\_trenddata**(datadf: *pandas.core.frame.DataFrame*, header: dict, param\_dico: dict) → dict

clinical main plots of a trend recordings

parameters df : pdDataframe

recorded data (MonitorTrend.data)

**header** [dict] recording parameters (MonitorTrend.header)

**param\_dico** [dict] plotting parameters (MonitorTrend.param)

**Returns** **afig\_dico**

**Return type** {names:fig\_obj}

**class** anesthplot.record\_main.**MonitorTrend**(filename: Optional[str] = None, load: bool = True)  
 Bases: anesthplot.record\_main.\_SlowWave

monitor trends recordings:

input = filename : path to file load = boolean to load data (default is True)

**file** [str] short name

**filename** [str] long name

**header** [dict] record parameters

**param** [dict] parameters

**clean\_trend** [external] clean the data

**show\_graphs** [external] plot clinical main plots

**class** anesthplot.record\_main.**TaphTrend**(filename: Optional[str] = None)  
 Bases: anesthplot.record\_main.\_SlowWave

taphonius trends recordings input : filename : path to file

data : pd.DataFrame = recorded data header : dictionary = recorded info (patient, ) param : dictionary  
 = usage information (file, scales, ) actions : pd.DataFrame

clean\_trend : to be developped show\_graphs : plot the clinical debrief suite

**extract\_events**()

**plot\_ventil\_drive**()

**plot\_events**(todrop: Optional[list] = None)

**export\_taph\_events**(save\_to\_file=False)  
 export in a txt files all the events (paths:~/temp/events.txt)

**class** anesthplot.record\_main.**TelevetWave**(filename=None)  
 Bases: anesthplot.record\_main.\_FastWave

class to organise teleVet recordings transformed to csv files. input:

filename : str (fullpath, default:None)

**class** anesthplot.record\_main.**MonitorWave**(filename: Optional[str] = None, load: bool = True)  
 Bases: anesthplot.record\_main.\_FastWave

**class to organise monitorWave recordings.** input : filename = path to file load = boolean to load data (default is True)

attibutes FILLME

methods FILLME

**anesthplot.record\_main.main**(file\_name: Optional[str] = None)

main script called from command line call : python anesthPlot/anesthplot/\_\_main\_\_.py args : optional filename (fullname)

return: set of plots for either monitorTrend, monitorWave oe televet recording



## MODULES

### 3.1 anesplot package

#### 3.1.1 Subpackages

##### anesplot.config package

##### Submodules

##### anesplot.config.build\_recordrc module

build a recordRc.yaml configuration file to adapt to a specific computer location at the root of anesplot

- input <-> data : to load the records
- output <-> save : to save the plots

---

`anesplot.config.build_recordrc.filedialog(kind="", directory="/Users/cdesbois/pg/chrisPg/anesthPlot/anesplot/config", for_open=True, fmt="", is_folder=False)`  
general dialog function.

`anesplot.config.build_recordrc.read_config()`  
locate & load the yaml file.

`anesplot.config.build_recordrc.write_configfile(path)`  
record the yaml file.

`anesplot.config.build_recordrc.main()`  
main function for script execution.

##### anesplot.config.load\_recordrc module

load an already generated recordRc.yaml configuration file

- input <-> data : to load the records
- output <-> save : to save the plots

`anesplot.config.load_recordrc.build_paths()`  
read the yaml configuration file.

`anesplot.config.load_recordrc.adapt_with_syspath(path_dico)`  
add the folder location to the system path.

## Module contents

### anesplot.loadrec package

#### Submodules

#### anesplot.loadrec.explore module

Created on Thu Mar 12 16:52:13 2020

@author: cdesbois

`anesplot.loadrec.explore.gui_choosefile(paths=None)`  
select a file via a dialog and return the file name.

#### anesplot.loadrec.loadmonitor\_trendrecord module

---

Created on Wed Jul 24 13:43:26 2019 @author: cdesbois

##### load a monitor trend recording:

- choose a file
- load the header to a dictionary
- load the data into a pandas dataframe

---

`anesplot.loadrec.loadmonitor_trendrecord.choosefile_gui(dirname: Optional[str] = None) → str`  
Select a file via a dialog and return the (full) filename.

**Parameters** `dir_path (str)` – location to place the gui (generally paths[data]) else home

**Returns** `fname[0]` – filename

**Return type** str

`anesplot.loadrec.loadmonitor_trendrecord.loadmonitor_trendheader(filename: str) → dict`  
load the file header.

**Parameters** `filename (str)` – full name of the file

**Returns** header

**Return type** dict

`anesplot.loadrec.loadmonitor_trendrecord.loadmonitor_trenddata(filename: str, headerdico: dict) → pandas.core.frame.DataFrame`  
load the monitor trend data

**Parameters**

- **filename** (*str*) – fullname
- **headerdico** (*dict*) – fileheader

**Returns** df = trends data

**Return type** pandas.DataFrame

### anesplot.loadrec.loadmonitor\_waverecord module

Created on Wed Jul 24 14:56:58 2019 @author: cdesbois

#### load a monitor wave recording:

- choose a file
- load the header to a pandas dataframe
- load the date into a pandas dataframe

anesplot.loadrec.loadmonitor\_waverecord.**choosefile\_gui**(*dirname: Optional[str] = None*) → str  
Select a file via a dialog and return the (full) filename.

**Parameters** **dir\_path** (*str*) – location to place the gui (generally paths[data]) else home

**Returns** **fname[0]** – filename

**Return type** str

anesplot.loadrec.loadmonitor\_waverecord.**loadmonitor\_waveheader**(*filename: Optional[str] = None*) → dict  
load the wave file header.

**Parameters** **filename** (*str*) – full name of the file

**Returns** header

**Return type** dictionary

anesplot.loadrec.loadmonitor\_waverecord.**loadmonitor\_wavedata**(*filename: Optional[str] = None*) → pandas.core.frame.DataFrame  
load the monitor wave csvDataFile.

**Parameters** **filename** (*str*) – full name of the file

**Returns** df = trends data

**Return type** pandas.DataFrame

### anesplot.loadrec.loadtaph\_trendrecord module

Created on Wed Jul 24 15:30:07 2019 @author: cdesbois

#### load a taphonius data recording:

- choose a file
- load the patient datafile to a dictionary
- load the physiological date into a pandas dataframe

**nb = 4 files per recording :**

- .pdf -> anesthesia record manual style
  - .xml -> taphonius technical record -> to be extracted
  - Patient.csv -> patient id and specifications
  - SDcsv -> anesthesia record
- 

`anesplot.loadrec.loadtaph_trendrecord.build_taph_decodedate_dico`(*pathdict: Optional[dict] = None*) → dict

list all the taph recordings and the paths to the record: input:

paths: dictionary containing {taph: pathToTheData}

**output:** dictionary: {date : filename}

`anesplot.loadrec.loadtaph_trendrecord.extract_record_day`(*monitor\_file\_name: str*) → str  
extract the date as YYYY\_MM\_DD from a monitor\_filename input:

monitor file name (shortname)

**output:** day : YYYY\_MM\_DD str

`anesplot.loadrec.loadtaph_trendrecord.choose_taph_record`(*monitorname: Optional[str] = None*) → str

select the taph recording: input:

taphdico : {date:path} builded from build\_taph\_decodedate\_dico() year = integer to place the pointer in pull down menu date = to be implemented (as year but to extract from monitor filename)

**output:** filename (str) full path

`anesplot.loadrec.loadtaph_trendrecord.loadtaph_trenddata`(*filename: str*) → pandas.core.frame.DataFrame

load the taphoniusData trends data.

**Parameters** `filename` (*str*) – fullname

**Returns** df = trends data

**Return type** pandas.DataFrame

`anesplot.loadrec.loadtaph_trendrecord.loadtaph_patientfile`(*filename: str*) → dict  
load the taphonius patient.csv file input:

**filename** [(*str*) the full filename] (the headername will be reconstructed inside the function)

**output:** descr = dict of patient\_data

## anesplot.loadrec.loadtelevet module

Created on Wed Jul 31 16:22:06 2019 @author: cdesbois

load televet exported (csv) data: to be developped

---

`anesplot.loadrec.loadtelevet.choosefile_gui(dirpath: Optional[str] = None) → str`  
select a file using a dialog.

**Parameters** `dir_path (str)` – optional location of the data (ex : paths[data], default : home)

**Returns** filename (full path)

**Return type** str

`anesplot.loadrec.loadtelevet.loadtelevet(fname: Optional[str] = None, all_traces: bool = False) → pandas.core.frame.DataFrame`

load the televetCsvExportedFile.

**Parameters**

- **file (str)** – name of the file
- **all\_traces (bool)** – load all the derivations

**Returns** df = recorded traces

**Return type** pandas.DataFrame

## Module contents

### anesplot.plot package

#### Submodules

### anesplot.plot.trend\_plot module

Created on Tue Apr 19 09:08:56 2016 @author: cdesbois

collection of functions to plot the trend data

---

`anesplot.plot.trend_plot.remove_outliers(df: pandas.core.frame.DataFrame, key: str, limits: Optional[dict] = None) → pandas.core.series.Series`

remove outliers input:

df : pandas.DataFrame key : a column label limits : dictionary of key:(limLow, limHigh)

**output:** pandas.series without the outliers

`anesplot.plot.trend_plot.color_axis(ax: matplotlib.axes._axes.Axes, spine: str = 'bottom', color: str = 'r')`  
change the color of the label & tick & spine.

**Parameters**

- **ax (matplotlib.pyplot.axis)** – the axis
- **spine (str)** – optional location in [bottom, left, top, right]

- **colors** (*str*) – optional color

```
anesplot.plot.trend_plot.append_loc_to_fig(ax: matplotlib.axes._axes.Axes, dt_list: list, label: str = 'g')
                                             → dict
```

append vertical lines to indicate a location eg: arterial blood gas

## Parameters

- **ax** (*matplotlib.pyplot.axis*) – the axis
- **dt\_list** (*[datetime]*) – list of datetime values
- **label** (*str*) – a key to add to the label (default is g)

`anesplot.plot.trend_plot.hist_co2_iso(data: pandas.core.frame.DataFrame, param: Optional[dict] = None) → matplotlib.figure.Figure`

CO2 and iso histogramme (NB CO2 should have been converted from % to mmHg)

**Parameters**

- **data** (*pandas.DataFrame*) – the trends recorded data
- **param** (*dict*) – dictionary of parameters

**Returns** `fig` *pyplot.figure*

`anesplot.plot.trend_plot.cardiovasc(data: pandas.core.frame.DataFrame, param: Optional[dict] = None) → matplotlib.figure.Figure`

cardiovascular plot

**Parameters**

- **data** (*pandas.DataFrame*) – the recorded trends data keys used :[ip1s, ip1m, ip1d, hr]
- **param** (*dict*) – dict(save: boolean, path[save], xmin, xmax, unit, dtime = boolean for time display in HH:MM format)

**Returns** `fig=` *pyplot.figure*

`anesplot.plot.trend_plot.cardiovasc_p1p2(data: pandas.core.frame.DataFrame, param: Optional[dict] = None) → pandas.core.frame.DataFrame`

cardiovascular plot with central venous pressure (p2)

**Parameters**

- **data** (*pandas.DataFrame*) – the trends recorded data keys used :[ip1s, ip1m, ip1d, hr, ip2s, ip2m, ip2d]
- **param** (*dict*) – dict(save: boolean, path[save], xmin, xmax, unit, dtime = boolean for time display in HH:MM format)

**Returns** `fig=` *pyplot.figure*

`anesplot.plot.trend_plot.co2iso(data: pandas.core.frame.DataFrame, param: Optional[dict] = None) → matplotlib.figure.Figure`

anesth plot (CO2/iso)

**Parameters**

- **data** (*pandas.DataFrame*) – the recorded data keys used :[ip1s, ip1m, ip1d, hr]
- **param** (*dictionary*) – dict(save: boolean, path[save], xmin, xmax, unit, dtime = boolean for time display in HH:MM format)

:returns `fig=` *pyplot.figure*

`anesplot.plot.trend_plot.func(ax, x, y1, y2, color='tab:blue', x0=38)`

`anesplot.plot.trend_plot.co2o2(data: pandas.core.frame.DataFrame, param: dict) → matplotlib.figure.Figure`

respiratory plot (CO2 and Iso)

**Parameters**

- **data** (*pandas.DataFrame*) – recorded trends data keys used :[ip1s, ip1m, ip1d, hr]
- **param** (*dict*) – dict(save: boolean, path[save], xmin, xmax, unit, dtime = boolean for time display in HH:MM format)

**Returns** `fig=` *pyplot.figure*

`anesplot.plot.trend_plot.ventil`(*data: pandas.core.frame.DataFrame, param=<class 'dict'>*) →  
 matplotlib.figure.Figure  
 plot ventilation parameters (.tvInsp, .pPeak, .pPlat, .peep, .minVexp, .co2RR, .co2exp )

#### Parameters

- **data** (*pandas.DataFrame*) – recorded data, keys used :[ip1s, ip1m, ip1d, hr]
- **param** (*dict*) – dict(save: boolean, path[save], xmin, xmax, unit, dtm = boolean for time display in HH:MM format)

**Returns** fig= pyplot.figure

`anesplot.plot.trend_plot.recruit`(*data: pandas.core.frame.DataFrame, param: dict*) →  
 matplotlib.figure.Figure  
 display a recruit manoeuver (.pPeak, .pPlat, .peep, .tvInsp)

#### Parameters

- **data** (*pandas.DataFrame*) – recorded data keys used :[ip1s, ip1m, ip1d, hr]
- **param** (*dict*) – dict(save: boolean, path[save], xmin, xmax, unit, dtm = boolean for time display in HH:MM format)

:returns fig= pyplot.figure

`anesplot.plot.trend_plot.ventil_cardio`(*data: pandas.core.frame.DataFrame, param: dict*) →  
 matplotlib.figure.Figure  
 build ventilation and cardiovascular plot

#### Parameters

- **data** (*pandas.DataFrame*) – teh recorded trends data keys used :[ip1s, ip1m, ip1d, hr]
- **param** (*dict*) – dict(save: boolean, path[save], xmin, xmax, unit, dtm = boolean for time display in HH:MM format)

**Returns** fig= pyplot.figure

`anesplot.plot.trend_plot.save_distri`(*data: pandas.core.frame.DataFrame, path: dict*)  
 save as **O\_.** the 4 distributions graphs for cardiovasc annd respi

`anesplot.plot.trend_plot.fig_memo`(*apath: str, fig\_name: str*)  
 append latex citation commands in a txt file inside the fig folder create the file iif it doesnt exist

## anesplot.plot.wave\_plot module

Created on Tue Apr 19 09:08:56 2016

@author: cdesbois

`anesplot.plot.wave_plot.color_axis`(*ax: matplotlib.axes.\_axes.Axes, spine: str = 'bottom', color: str = 'r'*)  
 change the color of the label & tick & spine.

#### Parameters

- **ax** (*matplotlib.pyplot.axis*) – the axis
- **spine** (*str*) – optional location in [bottom, left, top, right]
- **colors** (*str*) – optional color



`anesplot.plot.wave_plot.plot_wave(data: pandas.core.frame.DataFrame, keys: list, param: dict) → matplotlib.figure.Figure`

plot the waves recorded (from as5)

**Parameters**

- **data** (*pandas.DataFrame*) – the recorded trends data
- **keys** (*list*) – one or two in [wekg,ECG,wco2,wawp,wflow,wap]
- **{mini}** (*dict*) – limits in point value (index), maxi: limits in point value (index)}

**Returns fig** plt.figure the plot

**Returns lines** plt.line2D the line to animate

(Nb plot data/index, but the xscale is indicated as sec)

`anesplot.plot.wave_plot.get_roi(fig: matplotlib.figure.Figure, df: pandas.core.frame.DataFrame, params: dict) → dict`

use the drawn figure to extract the relevant data in order to build an animation

**Parameters waves** (*MonitorWave object*) – a wave recording

**Returns** a dictionary containing ylims, xlims(point, dtime and sec),

traces used to build the plot, the fig object :rtype: dictionary

**input:** fig : matplotlib.pyplot.figure df : dataframe used to build the figure params = paramter dictionary

**Returns** dictionary containing ylims, xlims(point, dtime and sec)

**Return type** roi

`anesplot.plot.wave_plot.create_video(data: pandas.core.frame.DataFrame, param: dict, roi: dict, speed: int = 1, save: bool = False, savename: str = 'example', savedir: str = '~')`

create a video from a figure input:

waves : waves object speed : integer, speed of the display save : boolean (default=False) savename : str (default=example) savedir : str (path, default=~

**Returns** .mp4 file .png file

## Module contents

Created on Tue Apr 19 09:08:56 2016

functions to plot the trend data

@author: cdesbois

## anesplot.treatrec package

### Submodules

#### anesplot.treatrec.clean\_data module

Created on Wed Jul 31 16:05:29 2019

@author: cdesbois

`anesplot.treatrec.clean_data.clean_trenddata(df)`  
remove artifacts in the recorded trends

#### anesplot.treatrec.ekg\_to\_hr module

Created on Wed Feb 12 16:52:00 2020 @author: cdesbois

function used to treat an EKG signal and extract the heart rate typically (copy, paste and execute line by line)

### 0. after

```
:: import pandas as pd
    import anesplot.record_main as rec from anesplot.treatrec import ekg_to_hr as tohr
```

### 1. load the data in a pandas dataframe:

(through classes `rec.MonitorTrend` & `rec.MonitorWave`)

```
trendname = '' # fullname
or
trendname = rec.choosefile_gui()
```

```
wavename = rec.trendname_to_wavename(trendname)
-
# load the data
trends = rec.MonitorTrend(trendname)
waves = rec.MonitorWave(wavename)
-
# format the name
name = trends.header['Patient Name'].title().replace(' ', '')
name = name[0].lower() + name[1:]
```

## 2. treat the ekg wave:

- get parameters
- build a dataframe to work with (waves)
- low pass filtering
- build the beat locations (beat based dataframe):

```
params = waves.param
ekg_df = pd.DataFrame(waves.data.wekg)
ekg_df['wekg_lowpass'] = rec.wf.fix_baseline_wander(ekg_df.wekg,
                                                    waves.param['fs'])
beat_df = tohr.detect_beats(ekg_df.wekg_lowpass, mult=1)
```

## 3. perform the manual adjustments required:

- based on a graphical display of beat locations, an rr values
- build a container for the manual corrections:

```
figure = tohr.plot_beats(ekg_df.wekg_lowpass, beat_df)
to_change_df = pd.DataFrame(columns=beat_df.columns.insert(0, 'action'))
```

- remove or add peaks : zoom on the figure to observe only one peak, then:

```
to_change_df = tohr.remove_beat(beat_df, ekg_df, to_change_df, figure)
or
to_change_df = tohr.append_beat(beat_df, ekg_df, to_change_df, figure,
                                yscale=1)
```

- combine to update the beat\_df with the manual changes:

```
beat_df = tohr.update_beat_df(beat_df, to_change_df,
                              path_to_file="", from_file=False)
```

- save the peaks locations:

```
tohr.save_beats(beat_df, to_change_df, savename='', savepath=None)
(# or reload
beat_df = pd.read_hdf('beatDf.hdf', key='beatDf') )
```

## 4. go from points values to continuous time:

```
beat_df = tohr.compute_rr(beat_df)
ahr_df = tohr.interpolate_rr(beat_df)
tohr.plot_rr(ahr_df, params)
```

## 5. append intantaneous heart rate to the initial data:

```
ekg_df = tohr.append_rr_and_ihr_to_wave(ekg_df, ahr_df)
waves.data = tohr.append_rr_and_ihr_to_wave(waves.data, ahr_df)
trends.data = tohr.append_ihr_to_trend(trends.data, waves.data, ekg_df)
```

## 6. save:

```
tohr.save_trends_data(trends.data, savename=name, savepath='data')
tohr.save_waves_data(waves.data, savename=name, savepath='data')
```

---

`anesplot.treatrec.ekg_to_hr.detect_beats(ser, fs=300, species='horse', mult=1)`  
detect the peak locations

### Parameters

- **ser** (*pandas.series*) – the data
- **fs** (*integer*) – sampling frequency
- **species** (*string*) – in [horse]
- **mult** (*float*) – correction / 1 for qRs amplitude

**Returns** df=pandas.DataFrame

`anesplot.treatrec.ekg_to_hr.plot_beats(ecg, beats)`  
plot ecg waveform + beat location

`anesplot.treatrec.ekg_to_hr.append_beat(beatdf, ekgdf, tochange_df, fig, lim=None, yscale=1)`  
locate the beat in the figure, append to a dataframe[toAppend]

### Parameters

- **beatdf** (*pandas.DataFrame*) – contains the point based location (p\_locs)
- **ekgdf** (*pandas dataframe*) – contains the wave recording ((wekg\_lowpass)
- **tochange\_df** (*pandas.DataFrame*) – to store the beats toAppend or toRemove
- **fig** (*pyplot.Figure*) – figure to find time limits
- **lim** (*integer*) – ptBasedLim optional to give it manually
- **yscale** (*float*) – amplitude mutliplication factor for detection (default=1)

**Returns** tochange\_df: incremented changedf (pt location)

**Return type** pandasDataFrame

methods :

locate the beat in the figure, append to a dataframe[toAppend] 0.: if not present : build a dataframe:

```
>>> to_change_df = pd.DataFrame(columns=['toAppend', 'toRemove'])
```

**1.: locate the extra beat in the figure (cf plot\_beats())** and zoom to observe only a negative peak

**2.: call the function:**

```
>>> to_change_df = remove_beat(beatdf, ekgdf, tochange_df, fig)
-> the beat parameters will be added the dataframe
```

**.in the end of the manual check, update the beat\_df**

- first : save beat\_df and to\_change\_df
- **second** [run:]

```
>>> beat_df = update_beat_df()
```

`anesplot.treatrec.ekg_to_hr.remove_beat(beatdf, ekgdf, tochange_df, fig, lim=None)`  
locate the beat in the figure, append to a dataframe[toRemove]

**0.: if not present build a dataframe:**

```
>>> to_change_df = pd.DataFrame(columns=['toAppend', 'toRemove'])
```

**1.: locate the extra beat in the figure (cf plot\_beats())** and zoom to observe only a negative peak

**2.: call the function:::**

```
>>> to_change_df = remove_beat(beatdf, ekgdf, tochange_df, fig)
-> the beat parameters will be added the dataframe
```

**.(in the end of the manual check, update the beat\_df**

- first : save beat\_df and to\_change\_df
- **second** [run]

```
>>> beat_df = update_beat_df()
```

`anesplot.treatrec.ekg_to_hr.save_beats(beatdf, tochangedf, savename="", dirpath=None)`  
save the beats locations as csv and hdf5 file

#### Parameters

- **beatde** (*pd.dataframes*) –
- **tochangedf** (*pandas.dataframe*) –
- **savename** (*filename*) –
- **dirpath** (*path to save in*) –
- **output** –
- ----- –
- **file** (*hdf*) –
- **key**='beatDf' –

`anesplot.treatrec.ekg_to_hr.update_beat_df(beatdf, tochangedf, path_to_file="", from_file=False)`  
implement in the beat location the manual corrections fromFile = True force the disk loading of the dataframes

`anesplot.treatrec.ekg_to_hr.compute_rr(beatdf, fs=None)`  
compute rr intervals (from pt to time)

#### Parameters

- **beatdf** (*pd.DataFrame*) – with p\_loc

- **fs** (*integer*) – sampling frequency

**Returns** with:  $rr = rr \text{ duration}$   $rrDiff = rrVariation$   $rrSqDiff = rrVariation^2$

**Return type** `pd.DataFrame`

`anesplot.treatrec.ekg_to_hr.interpolate_rr(beatdf, kind=None)`  
interpolate the beat\_df (pt -> time values)

**Parameters**

- **beatDf** (`pd.DataFrame`) –
- **kind** (*str*) – linear or cubic(default)

**Returns**  $espts$  = evenly spaced points  $rrInterpol$  = interpolated  $rr$

**Return type** `pd.DataFrame` with evenly spaced data

`anesplot.treatrec.ekg_to_hr.plot_rr(ahr_df, param, HR=False)`  
plot RR vs pt values +  $rrSqDiff$

**Parameters**

- **pdDataFrame** (*hr\_df* =) –
- **params** – dict containing  $fs$  as key

`anesplot.treatrec.ekg_to_hr.append_rr_and_ihr_to_wave(wave, ahrdf)`  
append  $rr$  and  $ihr$  to the waves based on  $pt$  value (ie index)

`anesplot.treatrec.ekg_to_hr.plot_agreement(trenddf)`  
plot  $ipIHR$  &  $ihr$  to check agreement

`anesplot.treatrec.ekg_to_hr.append_ihr_to_trend(trenddf, wavedf, ekgrdf)`  
append  $ihr$  (instantaneous heart rate) to the trends

`anesplot.treatrec.ekg_to_hr.save_trends_data(trenddf, savename="", dirpath='data')`  
save the trends data to a csv and hdf file, including an  $ihr$  column

**Parameters**

- **trenddf** (`pd.dataframes`) –
- **savename** (*str*) –
- **dirpath** (*str*) – path to save in (default= current working directory)
- **output** –
- -----
- **file** (*hdf*) –
- **key='trends\_data'** –

`anesplot.treatrec.ekg_to_hr.save_waves_data(wavedf, savename="", dirpath='data')`  
save the trends data to a hdf file, including an  $ihr$  column

**Parameters**

- **trenddf** (`pd.dataframes`) –
- **savename** (*str*) –  $dirpath$  : path to save in (default=data)
- **output** –
- -----

- **hdf\_file** –
- **key='waves\_data'** –

## anesplot.treatrec.extract\_hypotension module

Spyder Editor

This is a temporary script file.

`anesplot.treatrec.extract_hypotension.extract_hypotension(atrend, pamin=70)`  
 return a dataframe with the beginning and ending phses of hypotension

### Parameters

- **atrend** (*MonitorTrend object*) –
- **pamin** (*float= threshold de define hypotension on mean arterial pressure*) –
- **70** (*((default is)*) –

**Returns** **durdf** – transitions (up and down, in seconds from beginning) and duration in the hypotension state (in seconds)

**Return type** pandas DataFrame containing

`anesplot.treatrec.extract_hypotension.plot_hypotension(atrend, durdf, durmin=15, pamin=70)`  
 plot the hupotentions phases

### Parameters

- **atrend** (*TYPE*) – DESCRIPTION.
- **durdf** (*TYPE*) – DESCRIPTION.
- **durmin** (*TYPE, optional*) – DESCRIPTION. The default is 15.

**Returns** **fig** – DESCRIPTION.

**Return type** TYPE

`anesplot.treatrec.extract_hypotension.scatter_length_meanhypo(atrend, durdf)`  
 draw a scatter plot (hypotensive arterial value vs duration of hypotension) :param trends: :type trends: MonitorTrend :param durdf: :type durdf: pandas dataframe containing the value and duration

**Returns** **fig**

**Return type** matplotlib.pyplot figure

`anesplot.treatrec.extract_hypotension.plot_all_dir_hypo(dirname=None, scatter=False)`  
 walk throught the folder and plot the values

## anesplot.treatrec.hr\_to\_hrv module

`anesplot.treatrec.hr_to_hrv.build_hrv_limits(spec='horse')`  
return a dico containing HRV limits (VLF, LF, HF) input : spec in [horse, man]

## anesplot.treatrec.wave\_func module

Created on Fri Dec 8 12:46:41 2017

@author: cdesbois

`anesplot.treatrec.wave_func.fix_baseline_wander(data, fs=500)`  
BaselineWanderRemovalMedian.m from ecg-kit. Given a list of amplitude values (data) and sample rate (sr), it applies two median filters to data to compute the baseline. The returned result is the original data minus this computed baseline.

`anesplot.treatrec.wave_func.rol_mean(ser, win_lengh=1, fs=500)`  
returns a rolling mean of a RR serie

### Parameters

- **pd.Series** (*ser*) –
- **win\_lengh** (*integer*) – window lenght for averaging (in sec),
- **fs** (*int*) – sampling frequency

`anesplot.treatrec.wave_func.return_points(df, fig)`  
return a tuple containing the point values of ROI

### Parameters

- **df** (*anesthesia record dataframe*) –
- **fig** (*pyplot.figure*) –

### Returns ROI

**Return type** dict

`anesplot.treatrec.wave_func.restrict_time_area(dfl, mini=None, maxi=None)`  
return a new dataframe with reindexation

### Parameters

- **dfl** (*pandas.DataFrame*) –
- **mini** (*integer*) – miniPointValue
- **maxi** (*integer*) – maxiPointValue

**Return type** pandas.DataFrame



## Module contents

### 3.1.2 Submodules

### 3.1.3 Module contents

anesthPlot is a package to plot/use clinical anesthesia records for teaching

**three way to use it:**

1. **run directly anesplot from a terminal** -> PYTHONPATH=<pathToAnesthPlot> python -m anesplot -> generate a quick plotting of most interesting parts
2. **from an ipython terminal** -> import anesthPlot.anesplot.recordmain as rec -> trends = rec.MonitorTrend() -> waves = rec.MonitorWave() -> and use the objects trends and waves
3. import the module in a python environment (see below)

**(the presets are actually designed**

- for use with equine anesthesia
- to load data from a Monitor generated datex AS3/5 monitoring machine)

typical use when importing the module to build a clinical case

```
import os import sys
```

```
import numpy as np import pandas as pd
```

```
import anesplot.record_main as rec sys.path.append(os.path.expanduser("~/pg/utills")) from utills import saveGraph im-
port bloodGases.bgmain_manual as bgman
```

```
paths = rec.paths paths[save] = os.path.expanduser("~/toPlay/temp/") os.chdir(paths[save])
```

```
## globals def save_plot(name):
```

```
    filename = os.path.join(paths[save], fig, name) saveGraph(filename, ext=png, close=False, verbose=True)
```

```
def explore_hdf(filename):
```

```
    try: hdf = pd.HDFStore(filename) keys= [key.replace(/, ) for key in hdf.keys()] print( found h5_file { }
```

```
        that contains { }.format(filename, keys))
```

```
        hdf.close()
```

```
    except: print({ } is not an h5 file.format(filename))
```

```
saveName = os.path.join(paths[save], data, aname.h5)
```

```
explore_hdf(saveName)
```

```
## load and work trendName = rec.choosefile_gui(paths[data]) WaveName = rec.choosefile_gui(paths[data])
```

```
# build objects with headers trends = rec.MonitorTrend(trendName, load=True) waves = rec.MonitorWave(waveName,
load=True)
```

```
# or append data (pretreated ones) #trends.data = pd.read_hdf(saveName, trend_df) #waves.data =
pd.read_hdf(saveName, wave_df)
```

```
#remove filenames del waveName, trendName
```

```
# now you are ready to work with loaded trends and waves
```



## INDICES AND TABLES

- `genindex`
- `modindex`
- `search`



## PYTHON MODULE INDEX

### a

- `anesplot`, 19
- `anesplot.config`, 6
  - `anesplot.config.build_recordrc`, 5
  - `anesplot.config.load_recordrc`, 5
- `anesplot.loadrec`, 9
  - `anesplot.loadrec.explore`, 6
  - `anesplot.loadrec.loadmonitor_trendrecord`, 6
  - `anesplot.loadrec.loadmonitor_waverecord`, 7
  - `anesplot.loadrec.loadtaph_trendrecord`, 7
  - `anesplot.loadrec.loadtelevet`, 8
- `anesplot.plot`, 13
  - `anesplot.plot.trend_plot`, 9
  - `anesplot.plot.wave_plot`, 12
- `anesplot.record_main`, 3
- `anesplot.treatrec`, 19
  - `anesplot.treatrec.clean_data`, 13
  - `anesplot.treatrec.ekg_to_hr`, 13
  - `anesplot.treatrec.extract_hypotension`, 18
  - `anesplot.treatrec.hr_to_hrv`, 18
  - `anesplot.treatrec.wave_func`, 19



## A

[adapt\\_with\\_syspath\(\)](#) (in module *anesplot.config.load\_recordrc*), 6  
[anesplot](#)  
     module, 19  
[anesplot.config](#)  
     module, 6  
[anesplot.config.build\\_recordrc](#)  
     module, 5  
[anesplot.config.load\\_recordrc](#)  
     module, 5  
[anesplot.loadrec](#)  
     module, 9  
[anesplot.loadrec.explore](#)  
     module, 6  
[anesplot.loadrec.loadmonitor\\_trendrecord](#)  
     module, 6  
[anesplot.loadrec.loadmonitor\\_waverecord](#)  
     module, 7  
[anesplot.loadrec.loadtaph\\_trendrecord](#)  
     module, 7  
[anesplot.loadrec.loadtelevet](#)  
     module, 8  
[anesplot.plot](#)  
     module, 13  
[anesplot.plot.trend\\_plot](#)  
     module, 9  
[anesplot.plot.wave\\_plot](#)  
     module, 12  
[anesplot.record\\_main](#)  
     module, 3  
[anesplot.treatrec](#)  
     module, 19  
[anesplot.treatrec.clean\\_data](#)  
     module, 13  
[anesplot.treatrec.ekg\\_to\\_hr](#)  
     module, 13  
[anesplot.treatrec.extract\\_hypotension](#)  
     module, 18  
[anesplot.treatrec.hr\\_to\\_hrv](#)  
     module, 18  
[anesplot.treatrec.wave\\_func](#)

module, 19

[append\\_beat\(\)](#) (in module *anesplot.treatrec.ekg\_to\_hr*), 15  
[append\\_ihr\\_to\\_trend\(\)](#) (in module *anesplot.treatrec.ekg\_to\_hr*), 17  
[append\\_loc\\_to\\_fig\(\)](#) (in module *anesplot.plot.trend\_plot*), 9  
[append\\_rr\\_and\\_ihr\\_to\\_wave\(\)](#) (in module *anesplot.treatrec.ekg\_to\_hr*), 17

## B

[build\\_hrv\\_limits\(\)](#) (in module *anesplot.treatrec.hr\_to\_hrv*), 18  
[build\\_paths\(\)](#) (in module *anesplot.config.load\_recordrc*), 5  
[build\\_taph\\_decodedate\\_dico\(\)](#) (in module *anesplot.loadrec.loadtaph\_trendrecord*), 8

## C

[cardiovasc\(\)](#) (in module *anesplot.plot.trend\_plot*), 10  
[cardiovasc\\_p1p2\(\)](#) (in module *anesplot.plot.trend\_plot*), 10  
[choose\\_taph\\_record\(\)](#) (in module *anesplot.loadrec.loadtaph\_trendrecord*), 8  
[choosefile\\_gui\(\)](#) (in module *anesplot.loadrec.loadmonitor\_trendrecord*), 6  
[choosefile\\_gui\(\)](#) (in module *anesplot.loadrec.loadmonitor\_waverecord*), 7  
[choosefile\\_gui\(\)](#) (in module *anesplot.loadrec.loadtelevet*), 8  
[choosefile\\_gui\(\)](#) (in module *anesplot.record\_main*), 3  
[clean\\_trenddata\(\)](#) (in module *anesplot.treatrec.clean\_data*), 13  
[co2iso\(\)](#) (in module *anesplot.plot.trend\_plot*), 11  
[co2o2\(\)](#) (in module *anesplot.plot.trend\_plot*), 11  
[color\\_axis\(\)](#) (in module *anesplot.plot.trend\_plot*), 9  
[color\\_axis\(\)](#) (in module *anesplot.plot.wave\_plot*), 12  
[compute\\_rr\(\)](#) (in module *anesplot.treatrec.ekg\_to\_hr*), 16  
[create\\_video\(\)](#) (in module *anesplot.plot.wave\_plot*), 12

## D

`detect_beats()` (in module *anesplot.treatrec.ekg\_to\_hr*), 15

## E

`extract_hypotension()` (in module *anesplot.treatrec.extract\_hypotension*), 18

`extract_record_day()` (in module *anesplot.loadrec.loadtaph\_trendrecord*), 8

`extract_taph_actions()` (*anesplot.record\_main.TaphTrend* method), 4

## F

`fig_memo()` (in module *anesplot.plot.trend\_plot*), 12

`filedialog()` (in module *anesplot.config.build\_recordrc*), 5

`fix_baseline_wander()` (in module *anesplot.treatrec.wave\_func*), 19

`func()` (in module *anesplot.plot.trend\_plot*), 11

## G

`get_roi()` (in module *anesplot.plot.wave\_plot*), 12

`gui_choosefile()` (in module *anesplot.loadrec.explore*), 6

## H

`hist_cardio()` (in module *anesplot.plot.trend\_plot*), 10

`hist_co2_iso()` (in module *anesplot.plot.trend\_plot*), 10

## I

`interpolate_rr()` (in module *anesplot.treatrec.ekg\_to\_hr*), 17

## L

`loadmonitor_trenddata()` (in module *anesplot.loadrec.loadmonitor\_trendrecord*), 6

`loadmonitor_trendheader()` (in module *anesplot.loadrec.loadmonitor\_trendrecord*), 6

`loadmonitor_wavedata()` (in module *anesplot.loadrec.loadmonitor\_waverecord*), 7

`loadmonitor_waveheader()` (in module *anesplot.loadrec.loadmonitor\_waverecord*), 7

`loadtaph_patientfile()` (in module *anesplot.loadrec.loadtaph\_trendrecord*), 8

`loadtaph_trenddata()` (in module *anesplot.loadrec.loadtaph\_trendrecord*), 8

`loadtelevet()` (in module *anesplot.loadrec.loadtelevet*), 8

## M

`main()` (in module *anesplot.config.build\_recordrc*), 5

`main()` (in module *anesplot.record\_main*), 4

## module

*anesplot*, 19

*anesplot.config*, 6

*anesplot.config.build\_recordrc*, 5

*anesplot.config.load\_recordrc*, 5

*anesplot.loadrec*, 9

*anesplot.loadrec.explore*, 6

*anesplot.loadrec.loadmonitor\_trendrecord*, 6

*anesplot.loadrec.loadmonitor\_waverecord*, 7

*anesplot.loadrec.loadtaph\_trendrecord*, 7

*anesplot.loadrec.loadtelevet*, 8

*anesplot.plot*, 13

*anesplot.plot.trend\_plot*, 9

*anesplot.plot.wave\_plot*, 12

*anesplot.record\_main*, 3

*anesplot.treatrec*, 19

*anesplot.treatrec.clean\_data*, 13

*anesplot.treatrec.ekg\_to\_hr*, 13

*anesplot.treatrec.extract\_hypotension*, 18

*anesplot.treatrec.hr\_to\_hrv*, 18

*anesplot.treatrec.wave\_func*, 19

*MonitorTrend* (class in *anesplot.record\_main*), 4

*MonitorWave* (class in *anesplot.record\_main*), 4

## P

`plot_agreement()` (in module *anesplot.treatrec.ekg\_to\_hr*), 17

`plot_all_dir_hypo()` (in module *anesplot.treatrec.extract\_hypotension*), 18

`plot_beats()` (in module *anesplot.treatrec.ekg\_to\_hr*), 15

`plot_header()` (in module *anesplot.plot.trend\_plot*), 10

`plot_hypotension()` (in module *anesplot.treatrec.extract\_hypotension*), 18

`plot_one_over_time()` (in module *anesplot.plot.trend\_plot*), 10

`plot_rr()` (in module *anesplot.treatrec.ekg\_to\_hr*), 17

`plot_trenddata()` (in module *anesplot.record\_main*), 3

`plot_wave()` (in module *anesplot.plot.wave\_plot*), 12

## R

`read_config()` (in module *anesplot.config.build\_recordrc*), 5

`recrut()` (in module *anesplot.plot.trend\_plot*), 11

`remove_beat()` (in module *anesplot.treatrec.ekg\_to\_hr*), 16

`remove_outliers()` (in module *anesplot.plot.trend\_plot*), 9

`restrict_time_area()` (in module *anesplot.treatrec.wave\_func*), 19

`return_points()` (in module *anesplot.treatrec.wave\_func*), 19



`rol_mean()` (in module *anesplot.treatrec.wave\_func*), 19

## S

`save_beats()` (in module *anesplot.treatrec.ekg\_to\_hr*), 16

`save_distri()` (in module *anesplot.plot.trend\_plot*), 12

`save_graph()` (in module *anesplot.plot.trend\_plot*), 9

`save_trends_data()` (in module *anesplot.treatrec.ekg\_to\_hr*), 17

`save_waves_data()` (in module *anesplot.treatrec.ekg\_to\_hr*), 17

`scatter_length_meanhypo()` (in module *anesplot.treatrec.extract\_hypotension*), 18

`select_type()` (in module *anesplot.record\_main*), 3

`select_wave()` (in module *anesplot.record\_main*), 3

## T

`TaphTrend` (class in *anesplot.record\_main*), 4

`TelevetWave` (class in *anesplot.record\_main*), 4

`trendname_to_wavename()` (in module *anesplot.record\_main*), 3

## U

`update_beat_df()` (in module *anesplot.treatrec.ekg\_to\_hr*), 16

## V

`ventil()` (in module *anesplot.plot.trend\_plot*), 11

`ventil_cardio()` (in module *anesplot.plot.trend\_plot*), 11

## W

`write_configfile()` (in module *anesplot.config.build\_recordrc*), 5